



topspin

Bruker BioSpin



• **AU Programming**

TopSpin 2.1
Version 2.1.2

NMR Spectroscopy

think forward

Copyright (C) by Bruker BiosSpin GmbH

All rights reserved. No part of this publication may be reproduced, stored in a retrieval system, or transmitted, in any form, or by any means without prior consent of the publisher. Product names used are trademarks or registered trademarks of their respective holders.

This document was written by

NMR

(C); Bruker BioSpin GmbH

printed in Federal Republic

of Germany 03075-2008

Document No: SM/AU2.1.2

Document Part No: /01

Contents

Chapter 1	Introduction	5
1.1	What is new in TOPSPIN 2.1.	5
1.2	What is new in TOPSPIN 2.0.	5
1.3	What are AU programs?	6
1.4	Other Manuals describing AU Programs/Macros	6
1.5	Quick reference to using AU programs	7
1.6	Installing and compiling AU programs	7
1.7	Executing AU programs	8
1.8	Viewing AU programs	8
1.9	About AU macros	9
1.10	About Bruker library functions	9
1.11	Creating your own AU programs	10
1.12	How an AU program is translated into C-code	16
1.13	Listing of all predefined C-statements	19
Chapter 2	Inventory of AU macros and Bruker library functions	23
2.1	Naming conventions	23
2.2	Macros for dataset handling	25
2.3	Macros prompting the user for input.	27
2.4	Macros handling Topspin parameters.	27
2.5	Acquisition macros	29
2.6	Macros handling the shim unit and the sample changer	29
2.7	Macros handling the temperature unit	31
2.8	Macros handling the MAS unit	31
2.9	1D processing macros	33
2.10	Peak picking, integration and miscellaneous macros.	35
2.11	Macros for algebraic operations on datasets	35
2.12	Deconvolution macros	37
2.13	2D processing macros	37
2.14	Macros reading and writing projections etc.	40
2.15	3D processing macros	41
2.16	Spectral Width calculation macros.	42
2.17	Plot editor related macros.	42
2.18	Macros converting datasets	44
2.19	Macros to execute other AU programs, Topspin macros or commands	44
2.20	Bruker library functions.	45
2.21	Macros for loop control.	46
2.22	Macros to return from an AU program	46

Chapter 3	General AU macros	49
Chapter 4	TOPSPIN Interface functions	57
Chapter 5	Macros changing the current AU dataset	59
Chapter 6	Macros copying datasets	73
Chapter 7	Macros handling rows/columns	77
Chapter 8	Macros converting datasets	89
Chapter 9	Macros handling TOPSPIN parameters	95
Chapter 10	Macros for Plot Editor/autoplot	109
Chapter 11	Macros prompting the user for input.	119
Chapter 12	Bruker library functions	125
Chapter 13	List of Bruker AU programs	157
Chapter 14	TOPSPIN parameter types	175
	14.1 Integer parameters	176
	14.2 Float parameters	177
	14.3 Double parameters	178
	14.4 Character-string parameters	178

Chapter 1

Introduction

1.1 What is new in TOPSPIN 2.1

Changes in TOPSPIN 2.1 with respect to AU programs.

1. The macro XAU requires two arguments, allowing you to freely choose the arguments to be propagated from the calling AU program.
2. C-language argument syntax `iarg_v` and `iarg_c` can be used in AU programs.
3. New AU macros to delete data have been added: `DELETEPROCDATA`, `DELETEIMAGINARYDATA`, `DELETERAWDATA`, `DELETEPROCNO`, `DELETEEXPNO`, `DELETENAME`.
4. New AU macros to fetch/store nD data have been added: `FETCHPARN`, `FETCHPARNS`, `STOREPARN`, `STOREPARNS`.
5. The functions `getParfileDirForRead` and `getParfileDirForWrite` replace the functions `getstan` and `PathXWinNMR*`.

1.2 What is new in TOPSPIN 2.0

Changes in TOPSPIN 2.0 with respect to AU programs.

1. AU programs that contain a plotting command can be entered with the argument *noplot*. This argument prevents plotting.
2. All AU-macros, e.g. EF, APK, QUIT must be specified in capital letters. In previous versions of TOPSPIN and its predecessor XWIN-NMR, capital letters were recommended but not required.
3. New macros exist for automatic creation of Plot Editor layouts. Examples are LAYOUT_OBJ_1D and LAYOUT_ADD.

1.3 What are AU programs?

AU programs can be considered as user defined TOPSPIN commands. Any repetitive task is most effectively accomplished through an AU program. All commands which can be entered on the TOPSPIN command line can also be entered in an AU program in the form of macros. This includes selecting and changing datasets, reading and setting parameters, starting acquisitions, processing data and plotting the result. A simple AU program is nothing else than a sequence of such macros which execute the corresponding TOPSPIN commands. However, AU programs may also contain C-language statements. In fact, an AU program is a C-program because all AU macros are translated to C-statements. TOPSPIN automatically compiles AU programs to executable binaries, using a C-compiler.

TOPSPIN offers three other ways of creating user defined commands: TOPSPIN macros (not to be confused with AU macros), Tcl/Tk scripts and Python programs. They differ from AU programs in that they do not need to be compiled.

1.4 Other Manuals describing AU Programs/Macros

Creating and using AU programs is described and referred to in various other manuals:

- **Processing Reference Guide:** for each processing command for which an AU macro exists, this macro and its usage is specified.
- **Acquisition Reference Guide:** for each acquisition command for which an AU macro exists, this macro and its usage is specified.

- **NMR Guide:** AU programs can be sorted and listed according to their usage showing their names and short descriptions.
- **Plot Editor manual:** chapter about AU program macros for plotting.

1.5 Quick reference to using AU programs

Bruker delivers a library of standard AU programs with TOPSPIN. After TOPSPIN has been installed you must do the following in order to use them:

1. Run ***expinstall*** once to install all AU programs
2. Run ***compileall*** once to compile all AU programs
3. Enter the name of an AU program to execute it

Furthermore, you can write your own AU programs in the following way:

1. Enter ***edau <name>***
The file <name> will be opened with a text editor
2. Do one of the following:
 - Write your own AU program from scratch
 - Read in an existing AU program and modify it according to your needs
3. Click ***Save, exit and compile*** ¹.
4. Enter the name of the AU program to execute it.

After you have installed a new version of TOPSPIN, you must run ***expinstall*** and ***compileall*** again to install and compile both Bruker's and your own AU programs.

1.6 Installing and compiling AU programs

When you have installed a new version of TOPSPIN, you must install the library AU programs once by executing the TOPSPIN command ***expinstall***.

-
1. If you are not using the internal editor, you have to compile the AU program in a separate step with the command ***cplbruk <name>***.

Your own AU programs which you created under a previous version of TOPSPIN are still available, they only need to be re-compiled.

An AU program is automatically compiled, the first time it is executed, i.e. when its name is entered on the command line.

To compile an AU program without executing it:

- enter ***cplbruk <name>***

or

- enter ***edau <name>*** and click *Exit and Compile*.

To compile all Bruker AU programs:

- enter ***compileall***

1.7 Executing AU programs

Once an AU program has been installed, there are 3 different ways to execute it:

1. Enter the name of the AU program. This will work if:
 - The AU program is already compiled
 - No TOPSPIN command or macro ¹ with the same name exist
2. Enter ***edau***

A list of available AU programs will appear. Click on the AU program you want to execute and click *Execute*.

1.8 Viewing AU programs

You can view existing AU programs in the following

- a) Enter ***edau***

A dialog box listing all AU programs is opened. From the *Options* menu, you can choose to display *Bruker defined*, *User defined* or *All AU*

1. Here we refer to a TOPSPIN macro created with ***edmac***

programs.

- b) Click on an AU program in the list

When you select a Bruker AU program, it is shown in *view* mode which means you cannot edit it. When you click on a user-defined AU program it is shown in *edit* mode which means you can change it.

3. Enter *listall_au*

A list and a short description of all library AU programs is stored in the file `listall` in the users home directory. Note that this list is also available in Chapter 13 of this manual.

1.9 About AU macros

We will use the word *macro* rather often throughout this manual referring to AU macros. This should not be confused with TOPSPIN macros which are files containing a sequence of TOPSPIN commands. TOPSPIN macros are created with **edmac** and executed with **xmac**. An AU macro, however, is a statement in an AU program which defines one or more TOPSPIN commands, library functions or C-language statements. In its simplest form, an AU macro defines one TOPSPIN command. For example the macros ZG and FT execute the TOPSPIN commands **zg** and **ft**, respectively. Other macros like FETCHPAR and IEXPNO do not define TOPSPIN commands, their function is only relevant in the context of an AU program. More complex macros may contain several TOPSPIN commands and/or C-statements. All macros in AU programs should be written in capital letters. They are automatically translated to the corresponding C-code when the AU program is compiled. AU macros are defined in the file:

```
/tshome/prog/include/aucmd.h
```

1.10 About Bruker library functions

Bruker library functions are C-functions which are contained in Bruker libraries. They offer several features which are also used in the TOPSPIN interface, for example the display of a list of datasets from which the user can select one dataset. If you use a Bruker library function in an AU program the corresponding library is automatically linked to the AU program

during compilation. The most important and versatile Bruker library functions are described in chapter 9.

1.11 Creating your own AU programs

1.11.1 Writing a simple AU program

Before you start writing an AU program, you might want to check if an AU program already exists which (almost) meets your requirements. If this is not the case, you can write your own AU program in the following way:

1. Enter ***edau <au-name>***

Your preferred TOPSPIN text editor will be opened ¹

2. Do one of the following:

- Insert an existing library AU program and modify it to your needs.
- Write a new AU program using the macros as described in this manual.

The last macro in an AU program should always be QUIT (or QUITMSG).

3. Click *Save, exit and compile* ².

1.11.2 Using variables

Since AU programs are C programs you can use C-language variables. Several variables are already predefined for usage in AU programs. In fact, we distinguish three different types of variables: predefined dedicated variables, predefined general variables and user defined variables.

1.11.2.1 Predefined dedicated variables

Predefined dedicated variables have the following properties:

- they do not need to be declared in an AU program

1. To change the TOPSPIN text editor enter ***set*** and click *Miscellaneous*

2. If you are not using the internal editor, you have to compile the AU program in a separate step with the command ***cplbruk <name>***

- their declaration is automatically added during compilation
- they are known to the AU main body and to possible subroutines
- they are set implicitly by certain macros, e.g. the variable `expno` is set by macros like `DATASET` and `IEXPNO`
- they should not be set explicitly, so do NOT use statements like:

```
expno = 11;
FETCHPAR("NS", &expno)
```
- they can be evaluated in macros or C-statements, e.g.:

```
DATASET(name, expno, 2, disk, "guest")
i1=expno+1;
```
- examples of different types of predefined dedicated variables are:
 - char-string: `name`, `disk`, `user`, `name2`
 - integer: `expno`, `procno`, `loopcount1`, `loopcount2`, `lastparflag`

A complete list of all predefined dedicated variables with their types can be found in Chapter 1.13.2

1.11.2.2 Predefined general variables

Predefined general variables have the following properties:

- they do not need to be declared in an AU program
- their declaration is automatically added during compilation
- they are known to the AU main body but not to possible subroutines
- they can be freely used for various purposes
- examples of different types of predefined general variables are:

```
integer: i1, i2, i3
float: f1, f2, f3
double:d1, d2, d3
char-string: text
```

A complete list of all predefined general variables with their types and initial values can be found in Chapter 1.13.3.

1.11.2.3 User defined variables

For simple AU programs the number of predefined general variables is sufficient, you do not need to declare any additional variables. For more complex AU programs you might need more variables or you might want to use specific names. In these cases you can define your own variables in the AU program. User defined variables have the following properties:

- they must be declared at the beginning of an AU program
- they can be freely used for various purposes
- they are known to the main AU program but not to possible subroutines
- examples of declarations are:

```
int ivar1, ivar2;
float fvar1, fvar2, fvar3;
double dvar1, dvar2, dvar3;
char cstr1[20], cstr2[200];
```

1.11.3 Using AU macros with arguments

Several AU macros take one or more arguments. Arguments can be constants (values) or variables. In fact, an argument can be specified in four different ways as described here for the macro REXPNO:

- as a constant, e.g.:

```
REXPNO(3)
```

- as a predefined dedicated variable e.g.:

```
REXPNO(expno+1)
```

- as predefined general variable, e.g.:

```
i1=6;
REXPNO(i1)
```

- as a user defined variable, e.g.:

```
int my_exp;
....
my_exp=1;
REXPNO(my_exp)
```

It is very important that the arguments are of the correct type. Macros can take arguments of the type integer (like REXPNO), float, double or character-string.

Some macros, for example STOREPAR, take TOPSPIN parameters as arguments and each parameter is of a certain type. For example, the AU statement

```
STOREPAR("O1", d1)
```

stores the value of the variable d1 into the parameter O1. The predefined (double) variable d1 is used since O1 is of the type double. The second argument could also be a constant, e.g.:

```
STOREPAR("O1", 287.15)
```

A list of all TOPSPIN parameters and their type can be found in Chapter 14.

1.11.4 Using AU programs with arguments

An AU program can be used with arguments. Arguments are available within the AU program as C-languages variables:

`i_argc`: the number of arguments

`i_argv`: the arguments

`cmd`: all specified arguments concatenated

The first argument is the AU program pathname and the second argument always `exec`. So for an AU program entered without arguments, `i_argc = 2` and `cmd` is an empty string. For the example ***myau a1 a2***

```
i_argc = 4
i_argv[0] = myau
i_argv[1] = exec
i_argv[2] = a1
i_argv[3] = a2
cmd = "a1 a2"
```

Note that `cmd` is actually legacy code whose usage is discouraged. It may no longer be supported in future versions.

1.11.5 Using C-language statements

AU programs can contain AU macros but also C-language statements like:

- define statements, e.g.: #define MAXSIZE 32768
- include statements, e.g.: #include <time.h>
- variable declarations, e.g. int ivar;
- variable assignments, e.g.: ivar = 20;
- loop structures, e.g.: for, while, do
- control structures, e.g.: if-else
- C-functions, e.g.: strcpy, strcmp, sprintf

Important: several C-language statements (including declarations of variables) are already predefined and automatically added during compilation of the AU program.

An example of an AU program using macros and C-statements is:

```
int eno, pno;
char datapath [500], dataname[50], datauser[50], datadisk[200];
(void) strcpy (dataname,name);
(void) strcpy (datauser,user);
(void) strcpy (datadisk,disk);
eno = expno;
pno = procno;
(void) sprintf (datapath,"%s/data/%s/nmr/%s/%d/pdata/%d/title",
datadisk, datauser, dataname, eno, pno);
if ( (i1 = showfile (datapath)) < 0 )
{
    Proc_err (DEF_ERR_OPT,"Problems with showfile function");
}
QUIT
```

Note that QUIT is an AU macro, *strcpy* and *sprintf* are C-functions and *showfile* and *Proc_err* are Bruker library functions.

For an explanation of C-functions and more information on C-language we refer to the literature on C-programming.

1.11.6 Additional hints on C-statements

If you are using C-language code in your AU programs, then there are a few things to be considered.

1. Using C-language header files

Several C-language header files are automatically included in your AU program during compilation. If you are using C-code which requires additional header files you must write your AU program in a special way. The main AU program should be a call to a subroutine which performs the actual task of the AU program. The *include* statements for the header file must be entered between the main AU program and the subroutine. This gives the following structure:

```
subroutine(curdat, cmd)
QUIT
#include <headerfile.h>
subroutine(curdat,cmd)
char *curdat, *cmd;
{
    MACRO1
    MACRO2
}
```

Such a structure is used in several Bruker library AU programs (e.g. `amplstab`, `decon_t1`, etc.). Several Bruker library functions like `PrintExpTime`, `gethighest`, `getxwinvers`, `pow_next` and `unlinkpr` also require an include statement in the AU program (see Chapter 12).

2. Some macros, e.g. `IEXPNO` and `IPROCNO` change the current AU dataset but do not make it available for subsequent commands. If they are followed by a `CPR_exec` or any C-statement which access the current AU dataset, then you must precede that statement with `SETCURDATA` (see also the descriptions of `SETCURDATA`, `IEXPNO` etc. in Chapter 5).
3. If you are using C-languages loop statements like *for*, *do* or *while* or control statements like *if*, we strongly recommend to always put the body of such statements between `{}`. If the body only contains simple macros like `ZG` or `FT` you can omit them because these macro definitions already contain `{}`. However, more complex macros might internally

define C-statements that include loop or control structures. If such a macro is used within a loop or control structure in the AU program, then you create nested loops which require the usage of `{}`.

1.11.7 Viewing Bruker standard AU programs for macro syntax

The syntax of many AU macros is trivial, just enter the TOPSPIN command in capital letters. Other macros and especially Bruker library functions are more complex. A detailed description of frequently used AU macros and functions can be found in subsequent chapters of this manual. Alternatively, you can also look for an existing AU program containing this macro or function. If, for example, you want to know the syntax of the macro WRPA, just search for an AU program containing the text WRPA in the directory:

```
<ts_home>\prog\au\src.exam
```

using the Windows of Linux Search function.

1.12 How an AU program is translated into C-code

This paragraph is intended for users who want to get a deeper understanding of the compilation process. If you simply want to write and use AU programs you can skip this paragraph.

TOPSPIN automatically translates your AU program into C-language and compiles it. Files and directories used during AU program compilation are:

```
/<tshome>/exp/stan/nmr/au/makeau  
/<tshome>/exp/stan/nmr/au/vorspann  
/<tshome>/prog/include/aucmd.h  
/<tshome>/prog/include/inc
```

The compilation process is entirely controlled by the script `makeau` which performs the following steps.

1. The file `vorspann` is concatenated with your AU program. This file contains a variety of definitions including
 - the C-program *main* statement
 - *#include* statements of C-header files (which in turn contain other definitions)

- *#define* statements which define constants
 - predefined dedicated variables, e.g.: *name, disk, user, expno, procno*
 - predefined general variables, e.g. : *text, i1, i2, i3, f1, f2, f3, d1, d2, d3*
2. After `vorspann` and your code have been concatenated, all macro definitions are replaced according to their definitions as described in the file `aucmd.h` and in the `inc` directory. In some cases, the name of the macro is the name of one of the files in `inc` directory and the entire content of the file represents that macro.
 3. Step 2. results in a C program source file which this file is compiled and an executable program is created. By default, the compilation is done with the GNU C-compiler `gcc` which is delivered with TOPSPIN. The linking process is done with the native linker which is part of the native C-compiler `cc`. All AU program's source files reside in:

```
/<tshome>/exp/stan/nmr/au/src
```

executables will be stored into:

```
/<tshome>/prog/au/bin
```

The following section shows the result of concatenating `vorspann` with the following AU program:

```
EFP
APK
SREF
QUIT
```

For better presentation, only a part of `vorspann` is shown. All variables declared in `vorspann` are listed in chapter 1.10.

```
#include <stdio.h>
#include <stdlib.h>

.....

main(argc,argv)
int argc;
char **argv;
{
char curdat[PATH_MAX];
char arglist[BUFSIZ];
```

```

int modret;
...
(void)getcurdat(1, curdat, disk, user, type, name, &expno, &procno);
...
modret = AU_program(curdat,arglist);
}
.....

AU_program(curdat,cmd)
char *curdat;
char *cmd;
{
int i1=0,i2=0,i3=0;
float f1=0,f2=0,f3=0,f998=0,f999=0;
double d1=0,d2=0,d3=0;
char text[BUFSIZ/2];

GETCURDATA
EFP
APK
SREF
QUIT

```

Note that the macro QUIT defines the closing C-language '}' statement.

1.12.1 Using the native gcc compiler

By default, AU programs are compiled with the Bruker delivered gcc compiler. If you want to use the native operating system compiler, you can do that as follows.

1. From the Windows Explorer or LINUX file manager

Open the following file with a text editor:

```
<tshome>/exp/stan/nmr/au/makeau
```

2. Search for the following line:

```
# $opt_native = 1;
```

and remove the # character at the beginning of the line

3. Save and close the file.
4. Start TOPSPIN and compile your AU programs

Now, under Windows, the Visual C++ compiler will be used. Note that this is not a part of the standard operating system. Under LINUX, the default system GCC compiler is used.

To activate the native compiler for the current TOPSPIN session only, enter the following command:

```
env set DEBUG_MAKEAU=-native
```

1.13 Listing of all predefined C-statements

1.13.1 Including header files

The following C-language header files are automatically included during compilation:

stdio.h, stdlib.h, unistd.h, string.h, errno.h, math.h, limits.h, fcntl.h

which reside in the following directory:

Under Windows: <tshome>/GNU/usr/include

Under LINUX: /usr/include

and

errop.h, brukdef.h, lib/uni.h, lib/libcb.h, lib/util.h, sample.h, aucmd.h

which reside in the directory:

/tshome/prog/include

Note that under LINUX, the packages *glibc-kernheaders* and *glibc-devel* must be installed to be able to compile AU programs (see Installation Guide Linux).

1.13.2 Predefined dedicated variables

The following list contains all predefined dedicated variables, their type and the AU macros by which they are set. Note that most variables are set or modified by several macros and only one or two are listed here.

type	variable	set by macros
int	lastparflag	USELASTPARS, USECURPARS
int	loopcount1	TIMES/END
int	loopcount2	TIMES2/END
int	loopcount3	TIMES3/END
int	loopcountinf	TIMESINFINITE
char	disk[256]	GETCURDATA
char	user[64]	GETCURDATA
char	type[16]	GETCURDATA
char	name[64]	GETCURDATA
int	expno	GETCURDATA, IEXPNO
int	procno	GETCURDATA, IPROCNO
char	disk2[256]	GETCURDATA2
char	user2[64]	GETCURDATA2
char	type2[16]	GETCURDATA2
char	name2[64]	GETCURDATA2
int	expno2	GETCURDATA2
int	procno2	GETCURDATA2
char	disk3[256]	GETCURDATA3
char	user3[64]	GETCURDATA3
char	type3[16]	GETCURDATA3
char	name3[64]	GETCURDATA3
int	expno3	GETCURDATA3
int	procno3	GETCURDATA3
char	namelist[10][64]	SETDATASET
char	dulist[10][256]	SETDATASET
char	userlist[10][64]	SETDATASET

Table 1.1

type	variable	set by macros
char	parsetlist[10][16]	RPARSETLIST
char	pulproglis[10][32]	RPULPROGLIST
int	expnolist[15]	SETDATASET
int	procnolist[15]	SETDATASET
int	loopcountlist[15]	RLOOPCOUNTLIST
float	vtlist[128]	RVTLIST
int	xloopcount	ILOOPCOUNTLIST
int	xpulprog	IPULPROGLIST
int	xparset	IPARSETLIST
int	xdataset	IDATASETLIST
int	xvt	IVTLIST
int	listcount1	TIMESLIST
FILE	*textfilepointer	
FILE	*debug	
char	longpath[PATH_MAX]	
char	Hilfs_string[BUF- SIZ/2]	

Table 1.1**1.13.3 Predefined general variables**

The following list contains all predefined general variables, their types and

initial values:

type	variable	initial value
int	i1	0
int	i2	0
int	i3	0
double	d1	0
double	d2	0
double	d3	0
float	f1	0
float	f2	0
float	f3	0
float	f998	0
float	f999	0
char	text[BUFSIZ/2]	

Table 1.2

1.14 What to do after changing an AU program?

After changing a parameter in the userspecific AU program TopSpin must be updated to store the changed information in the AU program. This can be done with the command Show_meta (**SM_SHOWP**). This command also offers different arguments as following:

- **SM_RAW** ---- update raw data
- **SM_RAWP** ---- update acquisition parameters
- **SM_PROC** ---- updat processed data
- **SM_PROCP** ---- update processing parameters
- **SM_ALL** ----- update data and parameters
- **SM_SHOWR** ---- switch to raw data

- **SM_SHOWP** ---- switch to processed data
- **SM_DEL** ---- removed data
- **SM_PEAK** ---- update peaks
- **SM_INT** ---- update integrals

Please note that changing the peaklist with a macro in the AU program does not require the command **SM_PEAK**. The changings are implemented automatically.

Chapter 2

Inventory of AU macros and Bruker library functions

2.1 Naming conventions

This chapter lists most AU macros and Bruker library functions that are available for AU programming. Simple macros with their short description are only mentioned in this chapter. More complex macros and AU functions are mentioned here and described more extensively in the following chapters. Table 2.1 explains the macro conventions used in this chapter.

Macro	Explanation
XXX	The macro can be typed "as is". There is no further explanation for the macro in this manual.
XXX(arg1,arg2)	The macro XXX takes two arguments. Because the macro is easy to use, there is no further description in this manual.
XXX *	Like XXX, but there is a detailed description in one of the following chapters.
XXX(....) *	The macro XXX takes one or more arguments and its usage is described in one of the following chapters.

Table 2.1 Macro conventions

Several AU macros that are described in this chapter require one or more arguments. These arguments can be constants or variables as described in Chapter 1.11.3. It is very important to use the correct type of argument in a macro call. The macros described in the tables of this chapter use the following arguments:

integer : *i1*, *i2*, *i3*, *eno*, *pno*

float : *f1*

double : *d1*

char-string: *text*, *cmd*, *file*, *flag*, *mac*, *parm*, *parset*, *prog*, *shim*, *typ*, *dsk*,
usr, *nam*

Note that the arguments *i1*, *i2*, *i3*, *f1*, *d1* and *text* have the same names as the corresponding predefined general variables. The predefined general variables are easy to use because they do not need to be declared. You can, however, use your own variables as macro arguments.

2.2 Macros for dataset handling

Macro	Description
GETCURDATA *	Get the foreground dataset
SETCURDATA *	Make the current AU dataset available for subsequent AU statements
GETDATASET *	Prompt the user to specify a new dataset
DATASET(...)*	Set the current AU dataset
DATASET2(...)*	Set the 2nd dataset (like the TOPSPIN command edc2)
DATASET3(...)*	Set the 3rd dataset (like edc2)
GETCURDATA2	Read the 2nd dataset (like edc2)
GETCURDATA3	Read the 3rd dataset (like edc2)
DEXPNO *	Decrease the experiment number by one
IEXPNO *	Increase the experiment number by one
REXPNO(eno) *	Set the experiment number to the value of eno
DPROCNO *	Decrease the processing number by one
IPROCNO *	Increase the processing number by one
RPROCNO(i1) *	Set the processing number to the value of i1
GDATASETLIST	Prompt the user to enter a dataset list filename and read its contents
GLIST	Prompt the user to enter the dataset list filename and read its contents. In addition to the GDATASETLIST macro, GLIST also expects a pulse program and a parameter set name in the dataset list file.
DDATASETLIST	Decrement to the previous entry in the dataset list
IDATASETLIST	Increment to the next entry in the dataset list
RDATASETLIST(i1)	Read the dataset at position i1 of the dataset list and make it the current AU dataset

Macro	Description
IFEODATASETLIST	Checks if the end of the dataset list is reached. The answer is true if there is no further entry.
SETDATASET	Set the current AU dataset to the one currently defined by the dataset list
DU(dsk)	Set the disk unit (top level data directory) to <i>dsk</i>
SETUSER(usr)	Set the user name to the user <i>usr</i>
RE(name)	Read the dataset name.
WRA(eno) *	Copy the raw data to the experiment number <i>eno</i>
WRP(pno) *	Copy the processed data to the processing number <i>pno</i>
WRPA(...) *	Copy the raw and processed data to the specified dataset
VIEWDATA *	Show the current AU program dataset in a new window or activate the window that contains this dataset.
VIEWDATA_SAMEWIN*	Show the current AU program dataset in the current window.
AUDITCOMMENTA (cmt)	Add a user comment to the acquisition audit trail (audita.txt)
AUDITCOMMENTP(cmt)	Add a user comment to the processing audit trail (auditp.txt)
GDCHECK	generate checksum, making the processing audit trail consistent
GDCHECKRAW	generate checksum, making the raw audit trail consistent
ACQUPATH(x)	Returns the path of the file <i>x</i> in the acquisition data directory (ACQU)
PROCPATH(x)	Returns the path of the file <i>x</i> in the processed data directory (PROCNO)
DELETEPROCDATA	Delete processed data.

Macro	Description
DELETEIMAGINARY- DATA	Delete imaginary processed data.
DELETERAWDATA	Delete raw data.
DELETEPROCNO	Delete processed data directory (PROCNO).
DELETEEXPNO	Delete raw data directory (EXPNO).
DELETENAME	Delete data directory (NAME).

2.3 Macros prompting the user for input

Macro	Description
GETDOUBLE(text,d1) *	Prompt the user to enter a double value
GETFLOAT(text,f1) *	Prompt the user to enter a float value
GETINT(text,i1) *	Prompt the user to enter an integer value
GET- STRING(text,nam) *	Prompt the user to enter a text string

2.4 Macros handling Topspin parameters

Macro	Description
GETPROSOL *	Copy the probehead and solvent dependent parameters to the corresponding acquisition parameters
FETCHPAR(par,&val) *	Get an acquisition or processing parameter
FETCHPAR1(par,&val)	Get an F1 dimension parameter (2D acquisition/processing)
FETCHPAR3(par,&val)	Get an F1 dimension parameter (3D acquisition/processing)
FETCHPARS(par,&val) *	Get a status parameter (acquisition and processing)
FETCHPAR1S(par,&val)	Get an F1 dimension status parameter (2D)
FETCHPAR3S(par,&val)	Get an F1 dimension status parameter (3D)
FETCHPARN(dir,par,&val)	Get a parameter from specified direction (nD)
FETCHPARNS(dir,par,&val)	Get a status parameter from specified direction (nD)
STOREPAR(par,val) *	Store an acquisition, processing or output parameter
STOREPAR1(par,val)	Store an F1 dimension parameter (2D)
STOREPAR3(par,val)	Store an F1 dimension parameter (3D)

Macro	Description
STOREPARS(par,val) *	Store a status parameter (acquisition and processing)
STOREPAR1S(par,val)	Store an F1 dimension status parameter (2D)
STOREPAR3S(par,val)	Store an F1 dimension status parameter (3D)
STOREPARN(dir,par,&val)	Store a parameter to specified direction (nD)
STOREPARNS(dir,par,&val)	Store a status parameter to specified direction (nD)
FETCHPARAM(par,&val)	Get a tomography measurement parameter
STOREPARAM(par,val)	Store a tomography measurement parameter
FETCHT1PAR(par,&val)	Get a T1 parameter
STORET1PAR(par,val)	Store a T1 parameter
FETCHDOSYPAR(par,&val)	Get a dosy (eddosy) parameter
STOREDOSYPAR(par,val)	Store a dosy (eddosy) parameter
RPAR(parset,typ) *	Read a parameter set to the current dataset
WPAR(parset,typ) *	Write the current dataset parameters to a parameter set
DELPAR(parset)	Delete the parameter set <i>parset</i>

PARSET is not used in any AU program

2.5 Acquisition macros

Macro	Description
ZG	Start acquisition; if raw data already exist, they are overwritten
GO	Continue the acquisition on already existing raw data by adding to them
II	Initialize acquisition interface
RGA	Automatic receiver gain adjustment
MAKE_ZERO_FID	Create an empty FID
DEG90	Determine 90° pulse automatically
GPULPROGLIST	Prompt the user to enter the name of a pulse program list file and read its contents
DPULPROGLIST	Decrement to the previous name in the pulse program list
IPULPROGLIST	Increment to the next name in the pulse program list
RPULPRO- GLIST(i1)	Read the pulse program name in position i1 of the pulse program list and write it to the acquisition parameters
SETPULPROG	Store the current pulse program name from the pulse program list
IFEOPULPRO- GLIST	Check if the end of the pulse program list is reached. The answer is true if there is no further entry.

2.6 Macros handling the shim unit and the sample changer

Macro	Description
AUTOGAIN	Optimize lock gain
AUTOPHASE	Optimize lock phase
AUTOSHIM_ON	Turn autoshim on

Macro	Description
AUTOSHIM_OFF	Turn autoshim off
EJ	Eject sample from the magnet
IJ	Insert sample into the magnet
LOCK_ON	Turn lock on
LOCK_OFF	Turn lock off
ROT	Turn rotation on (use value RO from acquisition parameters)
ROTOFF	Turn rotation off and wait until rotation was turned off
LOPO	Set the lock parameters (lock power, lock gain, loop filter, loop time and loop gain)
LFILTER(i1)	Set the loop filter to the value of i1
LG	Auto-adjust the lock gain
LGAIN(f1)	Set the loop gain to the value of f1
LO(f1)	Set the lock power to the value of f1
LTIME(f1)	Set the loop time to the value of f1
LOCK	Lock according to the parameters LOCNUC and SOLVENT using the lock parameters from the edlock table
RSH(file)	Read the shim values from the specified file
SETSH(shim,i1)	Set one shim to the value of i1
SWEEP_ON	Turn the lock-sweep on
SWEEP_OFF	Turn the lock-sweep off
WSH(file)	Write the shim values to the specified file
TUNE(file)	Start autoshimming with the specified tune file
TUNESX	Start autoshimming with the tune file defined by the currently defined probehead and solvent

2.7 Macros handling the temperature unit

Macro	Description
TESET	Set the temperature on the temperature unit to the value of the acquisition parameter TE.
TEGET	Get the temperature from the temperature unit and store it in the acquisition status parameter TE
TE2SET	Set the temperature on the second regulator of the temperature unit to the value of the acquisition parameter TE2.
TE2GET	Get the temperature from the second regulator of the temperature unit and store it in the acquisition status parameter TE2
TEREADY(i1,f1)	After the temperature is set, wait until it is accurate to f1 degrees for at least 10 sec., then wait i1 seconds for stabilization.
TE2READY(i1,f1)	After the second temperature is set, wait until it is accurate to f1 degrees for at least 10 sec., then wait i1 seconds for stabilization.
TEPAR(file)	Read a file with parameter settings for the temperature unit
GVTLIST	Prompt the user to enter the variable temperature list name and read its contents
RVTLIST	Read the contents of the variable temperature list file defined by the acquisition parameter VTLIST.
DVTLIST	Decrement to the previous value in the vtlist
IVTLIST	Increment to the next value in the vtlist
VT	Read and set the temperature according to the current value of the vtlist

2.8 Macros handling the MAS unit

Macro	Description
MASE	Eject sample from MAS unit
MASI	Insert sample into MAS unit
MASR	Set spinning rate according to the acquisition parameter MASR
MASRGET	Get spinning rate from the MAS unit and store it in the status acquisition parameters
MASG(i1)	Start spinning of sample in MAS with at the most i1 retries
MASH	Halt spinning of sample in MAS

2.9 1D processing macros

Macro	Description
ABS	Automatic baseline correction (creates intrng file !)
ABSD	Automatic baseline correction with DISNMR algorithm (creates intrng file !)
ABSF	Automatic baseline correction between limits ABSF1 and ABSF2
APK	Automatic phase correction
APK0	Zero order automatic phase correction
APK1	First order automatic phase correction
APKF	Automatic phase correction using the spectral region determined by ABSF2 and ABSF1 for the calculation of the phase values.
APK0F	Zero order automatic phase correction using the spectral region determined by ABSF2 and ABSF1 for the calculation of the phase values.
APKS	Automatic phase correction especially suitable for polymer spectra
BC	Baseline correction of FID (DC correction)

Macro	Description
BCM	User defined spectrum baseline correction
CONVDTA(eno)	Convert digitally filtered FID into analogue (conventional) form
EF	Exponential window multiplication + Fourier transform
EFP	Exponential window multiplication + Fourier transform + phase correction using the processing parameters PHC0 and PHC1
EM	Exponential window multiplication of FID
FMC	Fourier Transform + magnitude calculation
FP	Fourier Transform + phase correction using the processing parameters PHC0 and PHC1
FT	Fourier Transform
GENFID(eno)	Create FID from processed data
GF	Gaussian window multiplication + Fourier Transform
GFP	Gaussian window multiplication + Fourier Transform + phase correction using the processing parameters PHC0 and PHC1
GM	Gaussian window multiplication
HT	Hilbert Transform
IFT	Inverse Fourier Transform
MC	Magnitude calculation
PK	Phase correction using the processing parameters PHC0 and PHC1
PS	Power spectrum calculation
QSIN	Squared sine window multiplication
SAB	Spline baseline correction using base_info file
SINM	Sine window multiplication
SINO	Calculate signal to noise ratio
SREF	Automatic spectral referencing using 2Hlock parameters
TM	Trapezoidal window multiplication

Macro	Description
TRF	Processing of the raw data according to the currently defined processing parameters
TRFP	Processing of the processed data according to the currently defined processing parameters
UWM	user-defined window multiplication

Note that 1D processing macros which access raw data, execute the corresponding command with the option **same**. For example, FT executes the command **ft same**.

2.10 Peak picking, integration and miscellaneous macros

Macro	Description
PP	Peak picking according to currently set processing parameters
PPH	Like PP, but with a peak histogram along the listing
PPP	Like PP, but the output is written to the file peaklist in the current processing data directory (PROCNO)
PPJ	Like PP, but store peaks in JCAMP-DX format
LI	List integrals according to the currently defined intrng file. The macro ABS can be used to create an intrng file.
LIPP	List integrals and all peaks in the integral ranges
LIPPF	Like LIPP, but works always on the full spectrum
PP2D	Perform peak picking on a 2D dataset.
RMISC(typ,file)	Read a file from one of the following list types: base_info, baslpnts, intrng, peaklist or reg
WMISC(typ,file)	Write a base_info, baslpnts, intrng, peaklist or reg file to its lists directory

2.11 Macros for algebraic operations on datasets

Macro	Description
ADD	Add 2nd and 3rd dataset and put the result into the current dataset. The 3rd dataset is multiplied by DC.
ADDFID	Add two FIDs multiplying one of them with DC.
ADDC	Add the constant DC to the current dataset
AND	Put logical "and" of 2nd and 3rd dataset into the current dataset
DIV	Divide 2nd and 3rd dataset and put the result into the current dataset. The 3rd dataset is multiplied by DC.
DT	Calculate the first derivative of the dataset
FILT	Apply a software digital filter to the current dataset

Macro	Description
LS	Left shift spectrum or FID by NSP points
MUL	Multiply 2nd and 3rd dataset and put the result into the current dataset. The 3rd dataset is multiplied by DC.
MULC	Multiply the current dataset with DC
NM	Negate current spectrum
RS	Right shift spectrum or FID by NSP points
RV	Reverse the spectrum
ZF	Zero the spectrum (1r,1i)
ZP	Zero the first NZP points of the spectrum or FID

2.12 Deconvolution macros

Macro	Description
GDCON	Gaussian deconvolution of the peaks automatically picked according to the currently set processing parameters
LDCON	Lorentzian deconvolution of the peaks automatically picked according to the currently set processing parameters
MDCON	Mixed Gaussian/Lorentzian deconvolution of the peaks in the peaklist file. The peaklist file can be created with the <i>ppp</i> command and it can be modified using the <i>edmisc</i> command.

2.13 2D processing macros

Macro	Description
ABS1	Baseline correction in F1 dimension
ABS2	Baseline correction in F2 dimension
ABSD1	Baseline correction in F1 dimension using the DISNMR algorithm
ABSD2	Baseline correction in F2 dimension using the DISNMR algorithm
ABSOT1	Trapezoidal baseline correction in F1 dimension using a slightly different algorithm than abst1
ABSOT2	Trapezoidal baseline correction in F2 dimension using a slightly different algorithm than abst2
ABST1	Trapezoidal baseline correction in F1 dimension using the processing parameters ABSF1, ABSF2, SIGF1, SIGF2
ABST2	Trapezoidal baseline correction in F2 dimension using the processing parameters ABSF1, ABSF2, SIGF1, SIGF2
ADD2D	Add the processed data of the 2nd dataset to the current dataset.
ADDSER	Add the raw data of the 2nd dataset to the current dataset.
BCM1	Baseline correction of all columns using the coefficients that were obtained with a manual 1D baseline correction
BCM2	Baseline correction of all rows using the coefficients that were obtained with a manual 1D baseline correction

Macro	Description
INVSF	Interchange the frequencies of the two dimensions
LEVCALC	Calculate the levels for the contour representation of the 2D matrix
PTILT	Tilt the 2D matrix by an arbitrary angle
PTILT1	Tilt the 2D matrix along its central vertical line
REV1	Reverse the spectrum in F1 dimension
REV2	Reverse the spectrum in F2 dimension
SUB1	Subtract 1D spectrum in F1 dimension (no change in sign !)
SUB2	Subtract 1D spectrum in F2 dimension (no change in sign !)
SUB1D1	Subtract 1D spectrum in F1 dimension
SUB1D2	Subtract 1D spectrum in F2 dimension
SYM	Symmetrize COSY spectrum
SYMA	Symmetrize phase sensitive COSY spectrum
SYMJ	Symmetrize J-resolved spectrum
TILT	Tilt J-resolved spectrum by an internally calculated angle
XF1	Fourier transform in F1 dimension
XF1P	Phase correction in F1 dimension using the processing parameters PHC0 and PHC1
XF2	Fourier transform in F2 dimension
XF2P	Phase correction in F2 dimension using the processing parameters PHC0 and PHC1
XFB	Fourier transform in both dimensions
XFBP	Phase correction in both dimensions
XF1M	Magnitude calculation in F1 dimension
XF2M	Magnitude calculation in F2 dimension
XFBM	Magnitude calculation in both dimensions
XF1PS	Power spectrum in F1 dimension
XF2PS	Power spectrum in F2 dimension
XFBPS	Power spectrum in both dimensions

Macro	Description
XHT1	Hilbert Transform in F1 dimension
XHT2	Hilbert transform in F2 dimension
XIF1	Inverse Fourier transform in F1 dimension
XIF2	Inverse Fourier transform in F2 dimension
XTRF	2D processing according to processing parameter flags (starts always on the raw data !)
XTRF2	2D processing according to F2 processing parameter flags only (starts always on the raw data !)
XTRFP	2D Processing according to the processing parameter flags
XTRFP1	2D processing according to the F1 processing parameter flags only
XTRFP2	2D processing according to the F2 processing parameter flags only
ZERT1	Zero a region of each column (F1). The region is determined by ABSF1/ABSF2 (first column) and SIGF1/SIGF2 (last column)
ZERT2	Zero a region of each row (F1). The region is determined by ABSF1/ABSF2 (first row) and SIGF1/SIGF2 (last row)
GENSER(en o)	Create a 2D series file from the processed data

Note that 2D processing macros which access raw data, execute the corresponding command with the option **same**. For example, XFB executes the command **xfb same**.

2.14 Macros reading and writing projections etc.

Macro	Description
F1SUM(i1,i2,pno)	Read sum of columns from i1 to i2 into the 1D processing number pno
F2SUM(i1,i2,pno)	Read sum of rows from i1 to i2 into the 1D processing number pno
F1DISCO(i1,i2,i3,pno)	Read disco projection between i1 and i2 columns with reference row i3 into the 1D processing number pno
F2DISCO(i1,i2,i3,pno)	Read disco projection between i1 and i2 rows with reference column i3 into the 1D processing number pno
F1PROJN(i1,i2,pno)	Read partial negative projection between columns i1 and i2 into the 1D processing number pno
F1PROJP(i1,i2,pno)	Read partial positive projection between columns i1 and i2 into the 1D processing number pno
F2PROJN(i1,i2,pno)	Read partial negative projection between rows i1 and i2 into the 1D processing number pno
F2PROJP(i1,i2,pno)	Read partial positive projection between rows i1 and i2 into the 1D processing number pno
RHNP(pno)	Read horizontal (F2) negative projection into the 1D processing number pno
RHPP(pno)	Read horizontal (F2) positive projection into the 1D processing number pno
RSC(i1,pno) *	Read column i1 of 2D into the 1D processing number pno
RSR(i1,pno) *	Read row i1 of 2D into the 1D processing number pno

Macro	Description
RVNP(pno) *	Read vertical (F1) negative projection into the 1D processing number pno
RVPP(pno) *	Read vertical (F1) positive projection into the 1D processing number pno
RSER(i1,eno,pno) *	Read row i1 of 2D raw data into the eno and pno
RSER2D(direc, i1,eno) *	Read plane number i1 in direction direc of 3D raw data into the eno
WSC(i1,pno,eno,nam,usr,dsk) *	Write a column back into position i1 of a 2D dataset defined by pno, eno, nam, usr and dsk
WSR(i1,pno,eno,nam,usr,dsk) *	Write a row back into position i1 of a 2D dataset defined by pno, eno, nam, usr and dsk.
WSER(i1,nam,eno,pno,dsk,usr)*	Write an FID back into position i1 of a 2D raw data defined by eno, pno, nam, dsk and usr.
WSERP(i1,nam,eno,pno,dsk,usr)	Write a processed FID back into position i1 of a 2D raw data defined by eno, pno, nam, dsk and usr.

2.15 3D processing macros

Macro	Description
TF3(flag,dsk)	Fourier transform in F3 dimension. The flag can be "y" or "n" and determines whether the imaginary parts are stored or not. The processed are stored on disk unit dsk.
TF2(flag)	Fourier transform in F2 dimension (flag as in TF3)
TF1(flag)	Fourier transform in F1 dimension (flag as in TF3)
TF3P(flag)	Phase correction in F3 dimension (flag as in TF3)
TF2P(flag)	Phase correction in F2 dimension (flag as in TF3)
TF1P(flag)	Phase correction in F1 dimension (flag as in TF3)
TABS3	Automatic baseline correction in F3 dimension

Macro	Description
TABS2	Automatic baseline correction in F2 dimension
TABS1	Automatic baseline correction in F1 dimension
R12(i1,pno)	Read F1-F2 plane into a new procno
R13(i1,pno)	Read F1-F3 plane into a new procno
R23(i1,pno)	Read F2-F3 plane into a new procno

2.16 Spectral Width calculation macros

Macro	Description
GETLIM	Get frequency of leftmost and rightmost peak from a 1D spectrum and adjust the sweep width of the 1D spectrum to the difference + 10%
GETLCOSY	Get frequency of leftmost and rightmost peak from a 1D spectrum and adjust the sweep width of a COSY spectrum to the difference + 10%
GETLXHCO	Get frequency of leftmost and rightmost peak from two 1D spectra and adjust the sweep width of an X-H correlation spectrum to the difference + 10%
GETLJRES	Get frequency of leftmost and rightmost peak from a 1D spectrum and adjust the sweep width of a J-RESolved spectrum to the difference + 10%
GETLINV	Get frequency of leftmost and rightmost peak from a 1D spectrum and adjust the sweep width of an INVerse spectrum to the difference + 10%

2.17 Plot editor related macros

XWP_LP *	Create a parameter listing for a plot with the plot editor.
XWP_PP *	Create a peak picking listing for a plot with the plot editor.

AUTO PLOT *	Plot the current dataset according to the plot editor layout defined by the processing parameter LAYOUT.
AUTO PLOT_TO_FILE(filename) *	as AUTO PLOT except that the plot is not sent to the printer but store in the postscript file <i>filename</i> .
DECLARE_PORTFOLIO *	Initialize the usage of other TOPSPIN portfolio AU macros. Required in XWIN-NMR and TOPSPIN ≤ 1.2 . Obsolete in TOPSPIN ≥ 1.3 .
CREATE_PORTFOLIO(filename) *	Create the TOPSPIN portfolio <i>filename</i> .
ADD_TO_PORTFOLIO(disk, user, name, expno, procno) *	Add the dataset that is specified with the arguments to the portfolio created with CREATE_PORTFOLIO.
ADD_CURDAT_TO_PORTFOLIO *	Add the current dataset to the portfolio created with CREATE_PORTFOLIO
CLOSE_PORTFOLIO *	Close the definition for the portfolio created with CREATE_PORTFOLIO. Must be used before AUTO PLOT_* macros.
AUTO PLOT_WITH_PORTFOLIO *	Plot the dataset(s) defined in the portfolio created with CREATE_PORTFOLIO according to the layout defined by the parameter LAYOUT.
AUTO PLOT_WITH_PORTFOLIO_TO_FILE(filename) *	as AUTO PLOT_WITH_PORTFOLIO except that the plot is not sent to the printer but store in the postscript file <i>filename</i> .
LAYOUT_ADD	Add object to current layout
LAYOUT_ADD_1D_OBJECT	Add 1D object to current layout
LAYOUT_ADD_PARAMETERS	Add parameter object to current layout
LAYOUT_BEGIN_FILE	Open layout file
LAYOUT_END_FILE	Close layout file

LAYOUT_END_FILE	Define layout format
LAYOUT_FORMAT	Define 1D object
LAYOUT_OBJ_1D	Define parameter object
LAYOUT_OBJ_PARAMETERS	Define title object
LAYOUT_OBJ_TITLE	Add object to current layout

2.18 Macros converting datasets

Macro	Description
FROMJDX(....) *	Convert a JCAMP-DX file to TOPSPIN data format
TOJDX(....) *	Convert a dataset to JCAMP-DX 6.0 format
TOJDX5(....) *	Convert a dataset to JCAMP-DX 5.0 format
JCONV(....) *	Convert a Jeol dataset to Bruker TOPSPIN format
VCONV(....) *	Convert a Varian dataset to Bruker TOPSPIN format

2.19 Macros to execute other AU programs, Topspin macros or commands

Macro	Description
CPR_exec(....) *	C-function for executing special TOPSPIN commands
WAIT_UNTIL(....) *	Hold the AU program until the specified date and time
XAUA	Execute the acquisition AU program stored in AUNM (<i>eda</i>). The next line in the AU program is executed after the AU program AUNM has finished.
XAUP	Execute the processing AU program stored in AUNMP (<i>edp</i>). The next line in the AU program is immediately executed after the AU program AUNMP has been started.

Macro	Description
XAUPW	Execute the processing AU program stored in AUNMP (edp). Like XAUP, but now the next line in the AU program is executed after the AU program AUNMP has finished.
XAU(prog, arg)	Execute the AU program prog with the wait option.
XCMD(cmd) *	Execute the TOPSPIN command for which no dedicated macro exists.
XMAC(mac)	Execute a TOPSPIN macro mac.

2.20 Bruker library functions

Macro	Description
CalcExpTime() *	Calculate the experiment time for the current experiment
PrintExpTime(...) *	Print the experiment time for the current experiment
check_pwd usr) *	Prompt the user usr to enter a password
GetNmrSuperUser() *	Get the name of the current TOPSPIN superuser
getdir(...) *	Get all file names and/or directory names within a directory
freedir(...) *	Free memory allocated by <code>getdir</code>
dircp(...) *	Copy a file
dircp_err(i1) *	Return the error message that corresponds to the error return value of a dircp function call
fetchstorpl(...) *	Read or store one or several plot parameters
FileSelect(...) *	Display a list of files and allow to select a file
gethighest(...) *	Return the next highest unused experiment number of a dataset
getParamDirs(...)	List all directories specified for key.
getParfileDirForRead(...) *	Determines pathname of list file to be read.

Macro	Description
getParfileDirForWrite(..) *	Determines pathname of list file to be written.
getstan(....) *	Return the pathname to the user's current experiment directory.
getxwinvers(....) *	Return the current version and patchlevel of TOPSPIN
mkudir(....) *	Create a complete directory path
PathXWinNMR() *	A class of functions which return pathnames to certain TOPSPIN directories
pow_next(i1) *	Round i1 to the next larger power of two
Proc_err(....) *	Show a message in a TOPSPIN dialog window
Show_status(text) *	Show a string in the status line of TOPSPIN
showfile(file) *	Show the contents of a file in a TOPSPIN window
ssleep(i1) *	Pause in an AU program for i1 seconds
unlinkpr(....) *	Delete all processed data files (1r, 1i, 2rr, 2ii etc.) of a dataset

2.21 Macros for loop control

Macro	Description
TIMES(n)	Execute the statements in the loop n times.
TIMES2(n)	Execute the statements in the loop n times. Normally used for the second level of nested loops.
TIMES3(n)	Execute the statements in the loop n times. Normally used for the second level of nested loops.
END	End of a loop.
STOP	Stop the AU program with the return value of AUERR.
STOPMSG("text")	Stop the AU program with the return value of AUERR and display the message "text"

2.22 Macros to return from an AU program

Macro	Description
ABORT	Abort the AU program or any of its subroutines with the return value of -1
ERRORABORT	Return from an AU program or any of its subroutines with the value of AUERR if it is less than 0
QUIT	Return from an AU program with the value of AUERR. QUIT is usually the last statement of the AU program code.
QUITMSG(text)	Print the text message and then return from the AU program with the value of AUERR. This is an alternative to QUIT.
STOP	Stop the AU program with the return value of AUERR.
STOPMSG("text")	Stop the AU program with the return value of AUERR and display the message "text"

Chapter 3

General AU macros

This chapter contains a description of all general AU macros which can be used for various purposes.

CPR_exec

NAME

CPR_exec - generic function for executing TOPSPIN commands

SYNTAX

CPR_exec(const char *command, int mode);

DESCRIPTION

CPR_exec is a generic function which can be used for executing TOPSPIN commands in AU. The first argument of CPR_exec is a string containing a Topspin command. The second argument must be one of the following values:

WAIT_TERM - wait for the command to finish, then start the next command

WAIT_START - wait for the command to start, then start the next command

CONT_EX - start the command and immediately start the next command

Practically all dedicated macros which execute a TOPSPIN command call CPR_exec with WAIT_TERM. For example, the macro FT is defined as

```
FT {SETCURDATA AUERR=CPR_exec("ft same",WAIT_TERM);}
```

The CPR_exec return value allows you to check for successful execution. The return value of CPR_Exec is NORM_TERM (=0) for normal termination or ERR_TERM (= -1) for error termination.

WAIT_START or CONT_EX can be used if asynchronous execution is required. For example, the AU macro XAUP uses WAIT_START to allow data simultaneous processing and acquisition in automation. Note that using WAIT_START and CONT_EX does not allow you to check the return value for successful execution.

For most commands a dedicated AU macro is available, like ZG for **zg** and FT for **ft**. If you want to use TOPSPIN commands for which no dedicated macro exist, e.g. editor commands or commands with special arguments, then you can use the generic macro XCMD which takes only one

argument, the TOPSPIN command and is started with WAIT_TERM.
XCMD is defined as:

```
XCMD(cmd) {SETCURDATA AU-  
ERR=CPR_exec(cmd,WAIT_TERM);}
```

In fact, the only reason to use CPR_exec explicitly is to start a command with WAIT_START or CONT_EX, i.e. to run commands simultaneously.

Note that dedicated macros and XCMD call SETCURDATA before they do their actual task. This ensures that they operate on the current AU dataset. If you use CPR_exec explicitly, it is recommended to precede it with SETCURDATA. Note that in the example below, CPR_exec is preceded by the macro ZG which implicitly calls SETCURDATA.

In summary:

- Use dedicated AU macros whenever you can
- Use XCMD when no dedicated macro is available
- Use CPR_EXEC when you want to use WAIT_START or CONT_EX

CPR_exec is part of the uni library which is delivered with TOPSPIN.

EXAMPLE

The following AU program gets the foreground dataset, runs an acquisition, starts the Fourier Transform and, after this has started, continues an acquisition on the next experiment number:

```
TIMES(10)  
ZG  
CPR_exec("ft", WAIT_START);  
IEXPNO  
END  
QUIT
```

SEE ALSO

XCMD - generic macro to execute commands for which no dedicated macro exists

SETCURDATA - make the current AU dataset available for subsequent AU statements

XAU

NAME

XAU - generic function for executing AU programs

SYNTAX

XAU(prog, arg)

DESCRIPTION

XAU is a general macro to execute (and, if necessary compile) AU programs. The macro takes two arguments:

prog - the AU program to be executed

arg - arguments

The second argument can be:

- any character string in "" containing one or more arguments
- cmd
to transfer all arguments from calling AU program
- ""
no arguments are propagated

In TOPSPIN 2.0 and older, XAU requires only one argument, the AU program to be executed, automatically propagates all arguments. In TOPSPIN 2.1 and newer, you can freely choose the arguments to be transferred.

User defined AU programs containing XAU macros must be modified. You could simply extend an XAU call with the extra argument cmd. As such, your AU program behaves exactly the same as in your previous version, namely propagating all caller arguments. The reason, however, that XAU was modified is that you normally do not want to propagate all arguments. In this case, you can replace XAU by XCMD, in which case the AU program can still be used with TOPSPIN 2.0 and older. In that case you can:

Replace XAU("auprog") by XCMD("auprog") or XAU("auprog", "")

to propagate no arguments

Replace XAU(auprog) by XAU("auprog", cmd)

to propagate all arguments

Use XCMD("auprog arg1 arg2") or XAU("auprog", "arg1 arg2")

to specify new arguments `arg1` and `arg2`.

SEE ALSO

XCMD - generic function for executing TOPSPIN commands

XCMD

NAME

XCMD - generic function for executing TOPSPIN commands

SYNTAX

XCMD(char *command)

DESCRIPTION

XCMD is a general macro to execute TOPSPIN commands for which no dedicated macro exists. For most TOPSPIN commands a dedicated macro does exist and we strongly recommend to:

Use dedicated macros whenever available

Note that XCMD executes CPR_exec with the option WAIT_TERM. If you want to use the options CONT_EX or WAIT_START, you must use CPR_exec.

If you want to check whether or not XCMD was executed successfully, you can check the value of AUERR (NORM_TERM or ERR_TERM).

EXAMPLE

The following AU program gets the foreground dataset, opens the acquisition parameter editor (**eda**) and runs an acquisition and Fourier transform:

```
XCMD("sendgui eda")
ZG
FT
QUIT
```

SEE ALSO

CPR_exec - C-function for executing special TOPSPIN commands.

WAIT_UNTIL

NAME

WAIT_UNTIL - hold the AU program until the specified date and time

SYNTAX

```
int WAIT_UNTIL(int hour, int minute, int day, int month)
```

DESCRIPTION

The function WAIT_UNTIL waits in an AU program until the specified date has been reached. The variables are internally converted to seconds. Every sixty seconds, the function checks whether the current date matches with the selected date. This function basically allows to program an event or command to start at a certain date rather than waiting for a certain time until something is executed.

EXAMPLE

Wait in the AU program until the 31st of October, 6 pm, and then continue:

```
WAIT_UNTIL(18,0,31,10)
```

SEE ALSO

ssleep - pause in an AU program for a certain number of seconds

Chapter 4

TOPSPIN Interface functions

AU programs are normally used to execute a series of acquisition or processing commands. For these commands you can use dedicated AU macros like ZG and FT. Less common is the use of TOPSPIN Java interface commands in AU programs. You can, however, do that with the XCMD or CPR_Exec macros and the command **sendgui**. Two examples are

```
XCMD("sendgui eda")  
to display the acquisition parameters
```

```
CPR_exec("sendgui .vr", WAIT_START);  
to perform a vertical reset of the current dataset
```

Note that XCMD is the same as CPR_exec with WAIT_TERM.

This can be used for all TOPSPIN interface commands like *Data window Tabs*,

Menu entries and Toolbar buttons. Here are some examples

TOPSPIN Interface	TOPSPIN Command	AU statement
<i>Menus</i>		
<i>File</i> ' <i>reopen</i>	reopen	XCMD("sendgui reopen")
<i>File</i> ' <i>Close</i>	close	XCMD("sendgui close")
<i>Window</i> ' <i>New Window</i>	newwin	XCMD("sendgui newwin")
<i>Data Window Tabs</i>		
<i>Spectrum</i>	spec	XCMD("sendgui spec")
<i>ProcPars</i>	edp	XCMD("sendgui edp")
<i>Title</i>	edti	XCMD("sendgui edti")
<i>Toolbar buttons</i>		
	.vr	XCMD("sendgui .vr")
	.zi	XCMD("sendgui .zi")
	.ov	XCMD("sendgui .ov")

Table 4.1

Chapter 5

Macros changing the current AU dataset

This chapter contains a description of all AU macros which can be used to change the current AU dataset, i.e. the dataset on which subsequent AU statements operate.

SETCURDATA

NAME

SETCURDATA - make the current AU dataset available for subsequent AU statements

SYNTAX

SETCURDATA

DESCRIPTION

SETCURDATA makes the current AU dataset, i.e. the dataset defined by the data path variables *disk*, *user*, *type*, *name*, *expno* and *procno*, available for subsequent AU commands. Normally, you do not need to enter SETCURDATA because it is automatically called by macros which operate on datasets before they perform their actual task. Furthermore, the macros DATASET and GETDATASET, which change the current AU dataset, automatically call SETCURDATA after they performed their actual task. In some cases, however, SETCURDATA must be specified explicitly in the AU program. For example, the macros IEXPNO and IPROCNO change the current AU dataset, but do not call SETCURDATA. If they are followed by a CPR_exec or any C-statement which access the current AU dataset, then you must precede that statement with SETCURDATA.

EXAMPLE

This example shows the part of the library AU program *multizg* which calculates the total experiment time of all acquisitions performed by this AU program:

```
int expTime;
static void PrintExpTime();
....
expTime = 0;
TIMES(i1)
    SETCURDATA;
    expTime += CalcExpTime() + 4;
```

```
IEXPNO;  
END  
DEXPNO;  
....  
QUIT
```

Note that IEXPNO is followed by SETCURDATA in the next cycle of the loop.

SEE ALSO

DATASET - set the current AU dataset

IEXPNO - increase the experiment number by one

DATASET

NAME

DATASET - set the current AU dataset

SYNTAX

DATASET(char *name, int expno, int procno, char *disk, char *user)

DESCRIPTION

The macro DATASET sets the current AU dataset. All data path variables *name*, *expno*, *procno*, *disk* and *user* must be specified as arguments. Subsequent AU commands will operate on this dataset.

EXAMPLE

The following AU program first gets the foreground dataset, then selects a new dataset and runs an acquisition:

```
char newname[20];
strcpy(newname, "glycerine");
DATASET(newname, expno, 3, disk, "peter")
ZG
QUIT
```

The data path variables in this example are entered in the following way:

- *expno* and *disk* keep the values of the current dataset
- *name* gets the value of *newname*, a variable defined in this AU program
- *procno* and *user* get the values *3* and *peter*, respectively, which are entered as constants

SEE ALSO

GETDATASET - prompt the user to specify a new dataset

DATASET2 - set the second dataset

IEXPNO - increase the experiment number by one

DATASET2, DATASET3

NAME

DATASET2 - set the second AU dataset
DATASET3 - set the third AU dataset

SYNTAX

DATASET2(char *name, int expno, int procno, char *disk, char *user)
DATASET3(char *name, int expno, int procno, char *disk, char *user)

DESCRIPTION

The macro DATASET2 sets the second AU dataset. The current (first) AU dataset is not affected by this macro. DATASET2 is typically used in combination with algebra macros, like ADD or MUL, which operate on the second and third dataset.

EXAMPLE

The following AU program gets the foreground dataset, adds the spectra of the next processing number and the one after that and stores the result into the current dataset:

```
DATASET2(name, expno, procno+1, disk, user)
DATASET3(name, expno, procno+2, disk, user)
ADD
QUIT
```

SEE ALSO

DATASET - set the current AU dataset
GETDATASET - prompt the user to specify a new dataset

GETDATASET

NAME

GETDATASET - prompt the user to specify a new dataset

SYNTAX

GETDATASET

DESCRIPTION

The macro GETDATASET prompts the user to specify a new dataset. A dialogue is opened and the user is requested to enter the data path variables *name*, *expno*, *procno*, *user* and *disk*. Subsequent AU commands will operate on this dataset. GETDATASET can be used anywhere in an AU program but, since it requires user input, should not be used in fully automated sequences.

NOTE

GETDATASET is not used very often. In AU programs, datasets are usually changed without user interaction, e.g. with the macros DATASET, IEXPNO etc.

EXAMPLE

The following AU program gets the foreground dataset, prompts the user to specify a new dataset and then processes this dataset:

```
GETDATASET
EFP
QUIT
```

SEE ALSO

DATASET - set the current AU dataset
IEXPNO - increase the experiment number by one
IPROCNO - increase the processing number by one

IEXPNO

NAME

IEXPNO - increase the experiment number by one

SYNTAX

IEXPNO

DESCRIPTION

The macro IEXPNO increases the experiment number of the current AU dataset by one. In fact, the value of the data path variable *expno* is incremented by one. Subsequent macros will operate on this new *expno*. IEXPNO is typically used in AU programs which run a series of acquisitions on datasets with the same *name* and successive *expnos*.

EXAMPLE

The following AU program gets the foreground dataset and runs acquisitions on eight successive *expnos*:

```
TIMES(8)
  ZG
  IEXPNO
END
QUIT
```

NOTE

IEXPNO must be followed by a SETCURDATA if the AU program continues with an explicit CPR_exec or C-statement (see SETCURDATA).

SEE ALSO

DEXPNO - decrease the experiment number by one
REXPNO - set the experiment number to the specified value
IPROCNO - increase the processing number by one
DATASET- set the current AU dataset

DEXPNO

NAME

DEXPNO - decrease the experiment number by one

SYNTAX

DEXPNO

DESCRIPTION

The macro DEXPNO decreases the experiment number of the current AU dataset by one. In fact, the value of the data path variable *expno* is decremented by one. Subsequent macros will operate on this new *expno*. DEXPNO is typically used after a loop which includes an IEXPNO at the end, to revert the effect of the last (unnecessary) IEXPNO.

EXAMPLE

The following AU program gets the foreground dataset, runs acquisitions on eight successive *expnos* and displays the data of the last *expno*:

```
TIMES(8)
  ZG
  IEXPNO
END
DEXPNO
VIEWDATA
QUIT
```

Note that DEXPNO must be followed by a SETCURDATA if the AU program continues with an explicit CPR_exec or C-statement (see SETCURDATA).

SEE ALSO

IEXPNO - increase the experiment number by one
REXPNO - set the experiment number to the specified value
DPROCNO - decrease the processing number by one

REXPNO

NAME

REXPNO - set the experiment number to the specified value

SYNTAX

REXPNO(int number)

DESCRIPTION

The macro REXPNO sets the experiment number of the current AU dataset to the specified value. In fact, the value of the data path variable *expno* is set. Subsequent macros will operate on this new *expno*.

EXAMPLE

The following AU program gets the foreground dataset, runs acquisitions on eight successive *expnos* then sets the current AU dataset back to the first *expno* and Fourier transforms it:

```
i1=expno;  
TIMES(8)  
  ZG  
  IEXPNO  
END  
REXPNO(i1)  
FT  
QUIT
```

Note that REXPNO must be followed by a SETCURDATA if the AU program continues with an explicit CPR_exec or C-statement (see SETCURDATA).

SEE ALSO

IEXPNO - increase the experiment number by one
DEXPNO - decrease the experiment number by one
RPROCNO - set the processing number to the specified value

IPROCNO

NAME

IPROCNO - increase the processing number by one

SYNTAX

IPROCNO

DESCRIPTION

The macro IPROCNO increases the processing number of the current AU dataset by one. In fact, the value of the data path variable *procno* is incremented by one. Subsequent macros will operate on this new *procno*. IPROCNO is typically used in an AU program which processes a series of datasets with same *name* and *expno* and successive *procnos*.

EXAMPLE

The following AU program runs Fourier transforms on eight successive *procnos*:

```
TIMES(8)
  FT
  IPROCNO
END
QUIT
```

Note that IPROCNO must be followed by a SETCURDATA if the AU program continues with an explicit CPR_exec or C-statement (see setcurdata).

SEE ALSO

DPROCNO - decrease the processing number by one

RPROCNO - set the processing number to the specified value

IEXPNO - increase the experiment number by one

DPROCNO

NAME

DPROCNO - decrease the processing number by one

SYNTAX

DPROCNO

DESCRIPTION

The macro DPROCNO decreases the processing number of the current AU dataset by one. In fact, the value of the data path variable *procno* is decremented by one. Subsequent macros will operate on this new *procno*. DPROCNO is typically used after a loop which includes an IPROCNO at the end, to revert the effect of the last (unnecessary) IPROCNO.

EXAMPLE

The following AU program gets the foreground dataset, runs a Fourier transform on eight successive *procnos* and displays the data of the last *procno*:

```
TIMES(8)
  FT
  IPROCNO
END
DPROCNO
VIEWDATA
QUIT
```

Note that DPROCNO must be followed by a SETCURDATA if the AU program continues with an explicit CPR_exec or C-statement (see setcurdata).

SEE ALSO

IPROCNO - decrease the experiment number by one
RPROCNO - set the processing number to specified value
DEXPNO - decrease the experiment number by one

RPROCNO

NAME

RPROCNO - set the processing number to the specified value

SYNTAX

RPROCNO(int number)

DESCRIPTION

The macro RPROCNO changes the current AU dataset by setting the processing number to the specified value. In fact, the value of the data path variable *procno* is set. Subsequent macros will then operate on this new *procno*.

EXAMPLE

The following AU program gets the foreground dataset and runs a Fourier transform on eight successive *procnos*. Then the current AU dataset is set back to the first *procno* which is then phase corrected:

```
i1=procno;  
TIMES(8)  
  FT  
  IPROCNO  
END  
RPROCNO(i1)  
APK  
QUIT
```

Note that RPROCNO must be followed by a SETCURDATA if the AU program continues with an explicit CPR_exec or C-statement (see SETCURDATA).

SEE ALSO

IPROCNO - increase the processing number by one
DPROCNO - decrease the processing number by one
REXPNO - set the experiment number to the specified value

VIEWDATA

NAME

VIEWDATA - Show the current AU program dataset in new window

SYNTAX

VIEWDATA

DESCRIPTION

The macro VIEWDATA shows the current AU program dataset in a new window or activates the window that contains this dataset. VIEWDATA is used whenever the current AU dataset is changed within the AU program, i.e. with DATASET, IEXPNO etc. and this dataset must be shown in TOPSPIN.

EXAMPLE

The following AU program gets the foreground dataset, increases the processing number and performs a Fourier transform storing the spectrum in this processing number. The spectrum is then shown in TOPSPIN:

```
IPROCNO  
FT  
VIEWDATA  
QUIT
```

SEE ALSO

VIEWDATA_SAMEWIN - show the current dataset in the current window
GETDATASET - prompt the user to specify a new dataset
DATASET - set the current AU dataset
IEXPNO - increase the experiment number by one
IPROCNO - increase the processing number by one

VIEWDATA_SAMEWIN

NAME

VIEWDATA_SAMEWIN - Show the current AU program dataset in the current window

SYNTAX

VIEWDATA_SAMEWIN

DESCRIPTION

The macro VIEWDATA_SAMEWIN shows the current AU program dataset in the current window that contains this dataset. It is used as an alternative to VIEWDATA.

EXAMPLE

The following AU program gets the foreground dataset, increases the processing number and performs a Fourier transform storing the spectrum in this processing number. The spectrum is then shown in TOPSPIN:

```
IPROCNO  
FT  
VIEWDATA_SAMEWIN  
QUIT
```

SEE ALSO

VIEWDATA - show the current dataset in new window
GETDATASET - prompt the user to specify a new dataset
DATASET - set the current AU dataset
IEXPNO - increase the experiment number by one
IPROCNO - increase the processing number by one

Chapter 6

Macros copying datasets

This chapter contains a description of all AU macros which can be used to copy the current AU dataset or parts of it to a new dataset.

WRA

NAME

WRA - copy the raw data to the specified experiment number

SYNTAX

WRA(int expno)

DESCRIPTION

The macro WRA copies the raw data, including the acquisition and processing parameters of the current AU dataset to a new experiment number. It does not copy the processed data.

EXAMPLE

The following AU program gets the foreground dataset and copies the raw data to eight successive experiment numbers, starting with *expno* 11:

```
i1=11;
TIMES(8)
  WRA(i1)
  i1++;
END
QUIT
```

SEE ALSO

WRP - copy the processed data to the specified processing number
WRPA - copy the raw and processed data to the specified dataset

WRP

NAME

WRP - copy the processed data to the specified processing number

SYNTAX

WRP(int procno)

DESCRIPTION

The macro WRP copies the processed data, including the processing parameters of the current AU dataset, to the specified processing number.

EXAMPLE

The following AU program gets the foreground dataset and copies the processed data to eight successive processing numbers, starting with *procno* 11:

```
i1=11;  
TIMES(8)  
    WRP(i1)  
    i1++;  
END  
QUIT
```

SEE ALSO

WRA - copy the raw data to the specified experiment number
WRPA - copy the raw and processed data to the specified dataset

WRPA

NAME

WRPA - copy the raw and processed data to the specified dataset

SYNTAX

WRPA(char *name, int expno, int procno, char *disk, char *user)

DESCRIPTION

The macro WRPA copies the raw and processed data of the current AU dataset to the specified dataset. WRPA takes 5 arguments, *name*, *expno*, *procno*, *disk* and *user*, i.e. the data path variables which define the dataset path. You can set one, several, or all of these variables to new values in order to define the destination dataset. You can, for instance, archive your data to an external medium by changing the value of the variable *disk* and leaving the other path variables the same.

EXAMPLE

The following AU program copies the current dataset to an external disk drive E:/:

```
WRPA(name, expno, procno, "E:/", user)
QUIT
```

SEE ALSO

WRA - copy the raw data to the specified experiment number
WRP - copy the processed data to the specified processing number

Chapter 7

Macros handling rows/columns

This chapter contains a description of all AU macros which can be used to read (write) rows or columns from (to) a 2D dataset and AU macros that can be used to read rows or planes from 3D raw data.

RSR

NAME

RSR - read a row from a 2D spectrum and store it as a 1D spectrum

SYNTAX

RSR(int row, int procno)

DESCRIPTION

The macro RSR reads a row from a 2D spectrum and stores it as a 1D spectrum. It can be used in the following ways:

- specified with procno > 0, executed on a 2D dataset
the specified row is stored under the current dataname, the current expno and the specified procno.
- specified with procno = -1, executed on a 2D dataset
the specified row is stored under dataset ~TEMP/1/pdata/1
- specified with procno > 0, executed on a 1D dataset
the specified row is read from a 2D dataset that resides under the current dataname, the current expno and the specified procno.
- specified with procno = -1, executed on a 1D dataset
the specified row is read from the 2D dataset from which the current 1D dataset was extracted (as defined in the file `used_from`).

EXAMPLE

The following AU program gets a 2D dataset and processes it. Then it reads row 16 and stores that under procno 999:

```
DATASET("my_2D_data", 1, 1, "C:/bio", "guest")
XFB
RSR(16, 999)
QUIT
```

SEE ALSO

RSC - read a column from a 2D spectrum and store it as a 1D spectrum

RSC

NAME

RSC - read a column from a 2D spectrum and store it as a 1D spectrum

SYNTAX

RSC(int column, int procno)

DESCRIPTION

The macro RSC reads a column from a 2D spectrum and stores it as a 1D spectrum. It can be used in the following ways:

- specified with procno > 0, executed on a 2D dataset
the specified column is stored under the current dataname, the current expno and the specified procno.
- specified with procno = -1, executed on a 2D dataset
the specified column is stored under dataset ~TEMP/1/pdata/1
- specified with procno > 0, executed on a 1D dataset
the specified column is read from a 2D dataset that resides under the current dataname, the current expno and the specified procno.
- specified with procno = -1, executed on a 1D dataset
the specified column is read from the 2D dataset from which the current 1D dataset was extracted (as defined in the file `used_from`).

EXAMPLE

The following AU program gets a 2D dataset and processes it in the F2 dimension. Then it reads column 128 and processes the resulting 1D dataset:

```
DATASET("my_2D_data", 1, 1, "C:/bio", "guest")
XF2
RSC(128, 10)
RPROCNO(10)
EF
QUIT
```

SEE ALSO

RSR - read a row from a 2D spectrum and store it as a 1D spectrum

WSC - replace a column of a 2D spectrum by a 1D spectrum

WSR

NAME

WSR - replace a row of a 2D spectrum by a 1D spectrum

SYNTAX

WSR(int row, int procno, int expno, char *name, char *user, char *disk)

DESCRIPTION

The macro WSR replaces a row of a 2D spectrum by a 1D spectrum. It can be used in the following ways:

- executed on a 1D dataset
the specified row of the specified dataset (must 2D data) is replaced by the current 1D data.
- executed on a 2D dataset
the specified row of the current 2D dataset is replaced by the specified dataset (must be 1D data)

EXAMPLE

The following AU program gets a 2D dataset, reads row 16, phase corrects this row and writes it back to the 2D data:

```
DATASET("my_2D_data", 1, 1, "C:/bio", "guest")
XFB
RSR(16, 999)
RPROCNO(999)
APK
WSR(16, 1, expno, name, user, disk)
QUIT
```

SEE ALSO

WSC - replace a column of a 2D spectrum by a 1D spectrum

RSR - read a row from a 2D spectrum and store it as a 1D spectrum

WSC

NAME

WSC - replace a column of a 2D spectrum by a 1D spectrum

SYNTAX

WSC(int column, int procno, int expno, char *name, char *user, char *disk)

DESCRIPTION

The macro WSC replaces a column of a 2D spectrum by a 1D spectrum. It can be used in the following ways:

- executed on a 1D dataset
the specified column of the specified dataset (must 2D data) is replaced by the current 1D data.
- executed on a 2D dataset
the specified column of the current 2D dataset is replaced by the specified dataset (must be 1D data)

EXAMPLE

The following AU program gets a 2D dataset, reads column 16, phase corrects this column and writes it back to the 2D data:

```
DATASET("my_2D_data", 1, 1, "C:/bio", "guest")
RSC(16, 999)
RPROCNO(999)
APK
WSC(16, 1, expno, name, user, disk)
QUIT
```

SEE ALSO

WSR - replace a row of a 2D spectrum by a 1D spectrum

RSC - read a column from a 2D spectrum and store it as a 1D spectrum

RSER

NAME

RSER - read a row from 2D or 3D raw data and store it as a 1D FID

SYNTAX

RSER(int row, int expno, int procno)

DESCRIPTION

The macro RSER reads a row from 2D or 3D raw data and stores it as a 1D fid. It can be used in the following ways:

- specified with expno > 0, executed on a 2D dataset
the specified row is stored under the current dataname and the specified expno. Processing parameters are stored under procno 1.
- specified with expno = -1, executed on a 2D dataset
the specified row is stored under dataset ~TEMP/1/pdata/1
- specified with expno > 0, executed on a 1D dataset
the specified row is read from a 2D raw data that resides under the current dataname and the specified expno. Processing parameters are read from procno 1.
- specified with expno = -1, executed on a 1D dataset
the specified row is read from the 2D dataset from which the current 1D dataset was extracted (as defined in the file `used_from`).

EXAMPLE

The following AU program splits 2D raw data into single fids that are stored in successive expnos:

```
int td;
FETCHPAR1S("TD",&td)
i1=0;
TIMES(td)
  i1 ++;
  RSER(i1,i1+expno,1)
```

```
END  
QUITMSG("--- splitter finished ---")
```

Note that this is the AU program ***splitser*** that is delivered with TOPSPIN.

SEE ALSO

WSER - replace a row of 2D raw data by 1D raw data

RSER2D - read a plane from 3D raw data and store it as 2D raw data

RSR - read a row from a 2D spectrum and store it as a 1D spectrum

WSER

NAME

WSER - replace a row of 2D raw data by 1D raw data

SYNTAX

WSER(int row, char *name, int expno, int procno, char *disk, char *user)

DESCRIPTION

The macro WSER replaces a row of 2D raw data by 1D raw data. It can be used in the following ways:

- executed on a 1D dataset
the specified row of the specified dataset (must be 2D data) is replaced by the current 1D data.
- executed on a 2D dataset
the specified row of the current 2D dataset is replaced by the specified dataset (must be 1D data)

EXAMPLE

The following AU program writes a number of 1D fids that are stored under the same data name and incremental expnos to 2D raw data.:

```
int ne, exp1, proc1;
char nm1[20];

ne=1; exp1=1; proc1=1;

strcpy(nm1, name);
GETSTRING("Enter name of 1D series: ", nm1)
GETINT("Enter starting EXPNO: ", exp1)
GETINT("Enter starting PROCNO: ", proc1)
GETINT("Enter # of Fids: ", ne)
USECURPARS
TIMES(ne)
    WSER(loopcount1+1, nm1, exp1, proc1, disk, user)
    exp1++;
END
```

QUIT

Note that this is the AU program *fidtoser* that is delivered with TOPSPIN.

SEE ALSO

RSER - read a row from 2D or 3D raw data and store it as a 1D FID

WSR - replace a row of a 2D spectrum by a 1D spectrum

WSC - replace a column of a 2D spectrum by a 1D spectrum

RSER2D

NAME

RSER2D - read a plane from 3D raw data and store it as 2D pseudo raw data

SYNTAX

RSER2D(char *direction, int plane, int expno, int procno)

DESCRIPTION

The macro RSER2D reads a plane from 3D raw data and stores it as 2D pseudo raw data. The first argument, the plane direction can be "s23" or "s13" for the F2-F3 or F1-F3 direction, respectively. The specified plane is stored under the current data name, the specified expno and the specified procno.

EXAMPLE

The following AU program gets a 3D dataset, reads the F2-F3-plane 64 and stores that under expno 11. It then switches to the output 2D dataset and processes it.

```
DATASET("my_3D_data", 1, 1, "C:/bio", "guest")
RSER2D("s23", 64, 11)
REXPNO(11)
XFB
QUIT
```

SEE ALSO

RSER - read a row from 2D or 3D raw data and store it as a 1D FID
WSER - replace a row of 2D raw data by 1D raw data

Chapter 8

Macros converting datasets

This chapter contains a description of all AU macros which can be used to convert TOPSPIN data. This includes the conversion of Bruker Aspect 2000/3000 data, Varian data and Jeol data to TOPSPIN data format as well as the conversion of TOPSPIN data to JCAMP-DX.

TOJDX, TOJDX5

NAME

TOJDX - convert a dataset to JCAMP-DX 6.0 format

TOJDX5 - convert a dataset to JCAMP-DX 5.0 format

SYNTAX

TOJDX(char *path, int type, int mode, char *title, char *origin, char *owner)

TOJDX5(char *path, int type, int mode, char *title, char *origin, char *owner)

DESCRIPTION

The macro TOJDX converts the current AU data to standard JCAMP-DX 6.0 format. It takes 6 arguments:

1. the pathname of the output file, e.g. /tmp/data1.dx
2. the output type: enter a number between 0 and 6, where:
 - 0 = FID (default)
 - 1 = real spectrum
 - 2 = complex spectrum
 - 3 = parameter files
 - 4 = raw data + real and imaginary processed data
 - 5 = raw data +real and imaginary processed data of all PROCNO's under the current EXPNO
 - 6 = raw data +real and imaginary processed data of all EXPNO's under the current NAME
3. the compression mode: enter 0, 1, 2 or 3
 - 0=FIX, 1=PACKED, 2=SQUEEZED, 3=DIFF/DUP (default)
4. the title as it appears in the output file: enter a character-string
5. the origin as it appears in the output file: enter a character-string
6. the owner as it appears in the output file: enter a character-string

If "*" is entered as an argument, then the default value is used.

Note that the macro TOJDX5 only supports the output types 0, 1, 2 and

3.

EXAMPLE

The following AU program gets the foreground dataset and performs a conversion to JCAMP on 5 successive experiment numbers. The name of the JCAMP file contains the *name* and *expno* of the corresponding TOPSPIN dataset.

```
TIMES(5)
  sprintf(text,"C:/TEMP/%s_%d.dx", name, expno);
  TOJDX(text, 0, 3, "*", "*", "*")
  IEXPNO
END
QUIT
```

SEE ALSO

FROMJDX - convert a JCAMP-DX file to TOPSPIN data format

FROMJDX

NAME

FROMJDX - convert a JCAMP-DX file to TOPSPIN data format

SYNTAX

```
FROMJDX(char *input-file)
```

DESCRIPTION

The macro FROMJDX converts a JCAMP-DX file to TOPSPIN data format. It takes one argument; the pathname of the input file, e.g. `/tmp/data1.dx`

FROMJDX can convert 1D and 2D data.

EXAMPLE

The following AU program converts all files with the extension `.dx` in the directory `C:/TEMP` to a TOPSPIN dataset:

```
char **listfile;
i1 = getdir ("C:/TEMP",&listfile,"*.dx");
TIMES(i1)
    sprintf(text, "C:/TEMP/%s", listfile[loopcount1]);
    FROMJDX(text)
END
QUIT
```

SEE ALSO

TOJDX - convert a dataset to JCAMP-DX format

`getdir` - get all file names and/or directory names within a directory

VCONV

NAME

VCONV - convert a Varian dataset to Bruker TOPSPIN format

SYNTAX

VCONV(char *v-name, char *x-name, int expno, char *disk, char *user)

DESCRIPTION

The macro VCONV converts a Varian dataset to TOPSPIN data format. It takes 5 parameters:

1. the name of the input Varian dataset
2. the name of the output TOPSPIN dataset
3. the experiment number of the output TOPSPIN dataset
4. the disk unit of the output TOPSPIN dataset
5. the user of the output TOPSPIN dataset

EXAMPLE

The following AU program converts a Varian dataset to TOPSPIN format:

```
VCONV("pinen_h.fid", "pinen_h", 1, "C:/bio", "joe")  
QUIT
```

Note that VCONV searches for the input data file in the directory defined by the environment variable VNMR. Assume the file resides in C:/bio. You can set VNMR from the TOPSPIN command line with:

```
env set VNMR=c:/bio
```

or inside the AU program with:

```
CPR_exec("env set VNMR=C:/bio", WAIT_TERM);
```

SEE ALSO

JCONV - convert a Jeol dataset to Bruker TOPSPIN format

JCONV

NAME

JCONV - convert a Jeol dataset to Bruker TOPSPIN format

SYNTAX

JCONV(char *j-name, char *x-name, int expno, char *disk, char *user)

DESCRIPTION

The macro JCONV converts a Jeol dataset to TOPSPIN data format. It takes 5 parameters:

1. the name of the input Jeol dataset
2. the name of the output TOPSPIN dataset
3. the experiment number of the output TOPSPIN dataset
4. the disk unit of the output TOPSPIN dataset
5. the user of the output TOPSPIN dataset

Note that JCONV searches for the input data file in the directory defined by the environment variable JNMR. Assume the file resides in C:/bio. You can set JNMR from the TOPSPIN command line with:

env set JNMR=c:/bio

or inside the AU program with:

```
CPR_exec("env set JNMR=C:/bio", WAIT_TERM);
```

EXAMPLE

The following AU program converts a Jeol dataset to TOPSPIN format:

```
JCONV("gx400h.gxd", "gx400h", 1, "C:/bio", "joe")  
QUIT
```

SEE ALSO

VCONV - convert a Varian dataset to Bruker TOPSPIN format

Chapter 9

Macros handling TOPSPIN parameters

This chapter contains a description of AU macros which can be used to get and store TOPSPIN parameters. Parameters are subdivided in acquisition, processing, output and plot parameters. Furthermore, they exist in two different forms; as foreground and status parameters. Finally, multi-dimensional datasets have parameter sets for each dimension. Different AU macros are available for getting and storing parameters of all categories, forms or dimensions.

FETCHPAR

NAME

FETCHPAR - get an acquisition, processing or output parameter

SYNTAX

FETCHPAR(par, &val)

DESCRIPTION

The macro FETCHPAR gets the value of a foreground parameter and stores it into an AU variable. This AU variable can then be used in subsequent AU statements. FETCHPAR allows to get acquisition parameters (*eda*) and processing parameters (*edp*). It is typically used to check or modify a parameter prior to an acquisition or processing statement.

The macro FETCHPAR takes two arguments:

1. the name of the parameter
2. the AU variable into which the parameter value will be stored

There are two important things to be considered:

1. The type of the AU variable must be the same as the type of the parameter (see Chapter 14).
2. The second argument must be specified as the variable's address, i.e. it must be prepended with the '&' character. This, however, does not count for a text variable since a text variable is already an address.

FETCHPAR works on 1D, 2D or 3D datasets and always gets a parameter of the first dimension (F2 for 1D, F2 for 2D and F3 for 3D).

The handling of the macros FETCHPAR1, FETCHPAR3, FETCHPARM, FETCHT1PAR and FETCHDOSYPAR is equivalent to the handling of FETCHPAR.

EXAMPLES

The following AU program gets the value of the processing parameter SI and processes the data 4 times, each time doubling the spectrum size

and storing the data in successive processing numbers:

```
FETCHPAR("SI", &i1)
TIMES(4)
  EFP
  IPROCNO
  i1 = i1*2;
  STOREPAR("SI", i1)
END
QUIT
```

The following AU statements get the values of the acquisition parameter DW and the processing parameter STSI and stores them in the predefined variables *f1* and *i1*, respectively. Then it gets value of the parameter ABSF1 and stores it in the user defined variable *leftlimit*.

```
float leftlimit;
...
FETCHPAR("DW", &f1)
FETCHPAR("STSI", &i1)
FETCHPAR("ABSF1", &leftlimit)
```

SEE ALSO

FETCHPARS - get a status parameter
FETCHPARN - get a parameter from specified direction
STOREPAR - store an acquisition, processing or output parameter

FETCHPARS

NAME

FETCHPARS - get a status parameter (acquisition and processing)

SYNTAX

FETCHPARS(par, &val)

DESCRIPTION

The macro `FETCHPARS` gets the value of a status parameter and stores it into an AU variable. This AU variable can then be used in subsequent AU statements. Acquisition status parameters are set by acquisition commands and describe the status of the dataset after acquisition. Note that the status parameters (*dpa*) describe what really happened and that this is sometimes different from what was set up before the acquisition as acquisition parameters (*eda*). For example, the status NS is smaller than originally specified when an acquisition was halted prematurely. Any AU program statement which follows an acquisition command and evaluates acquisition parameters must read status parameters. Therefore, `FETCHPARS` is typically used after acquisition or processing statements, for example for error or abort conditions (see example below).

The macro `FETCHPARS` takes two arguments:

1. the name of the parameter
2. the AU variable into which the value is value will be stored

There are two important things to be considered:

1. The type of the AU variable must be the same as the type of the parameter (see Chapter 14).
2. The second argument must be specified as the variable's address, i.e. it must be prepended with the '&' character. This, however, does not count for a text variable since a text variable is already an address.

The handling of the macros `FETCHPARS1` and `FETCHPARS3` is equivalent to the handling of `FETCHPARS`.

EXAMPLE

The following AU program performs a series of acquisitions on the same dataset until a minimum signal/noise is reached. In a loop 8 scans are acquired, Fourier transformed and phase corrected. Then the signal/noise of the spectrum is calculated and compared with the minimum value. If the minimum signal/noise was not reached yet, 8 more scans are accumulated etc. A maximum of 8000 scans is acquired. After the acquisition has been stopped, the total number of actually acquired scans is displayed.

```
STOREPAR("NS", 8)
GETFLOAT("Please enter the minimum signal/noise", f1)
ZG
TIMES(1000)
  FT
  APK
  SINO
  FETCHPARS("SINO", f2)
  if (f1 >= f2)
    break;
GO
END
FETCHPARS("NS", i1)
Proc_err (DEF_ERR_OPT,"Acquisition stopped after %d scans", i1);
QUIT
```

SEE ALSO

FETCHPAR - get an acquisition, processing or output parameter
FETCHPARNS - get a status parameter from specified direction
STOREPARS - store a status parameter (acquisition and processing)

STOREPAR

NAME

STOREPAR - store an acquisition, processing or output parameter

SYNTAX

STOREPAR(par, val)

DESCRIPTION

The macro STOREPAR stores the value of an AU variable into a parameter. This AU variable can then be used in subsequent AU statements. STOREPAR can be used for acquisition parameters (**eda**) and processing parameters (**edp**). It is typically used to set parameters prior to an acquisition or processing statement. STOREPAR takes two arguments:

1. the name of the parameter
2. the value to be stored which can be specified in two different forms:
 - as a constant
 - as the name of an AU variable

Important: the type of the parameter must be the same as the type of the constant or variable. (see Chapter 14).

NOTES

STOREPAR works on 1D, 2D or 3D datasets and always stores a parameter of the first dimension (F2 for 1D, F2 for 2D and F3 for 3D).

The handling of the macros STOREPAR1, STOREPAR3, STORET1PAR and STOREDOSYPAR is equivalent to the handling of STOREPAR.

EXAMPLE

The following AU program reads a standard parameter set, sets the pulse program and power level and asks the user for the number of scans. Then a dataset is acquired and processed according to these parameters.

```
RPAR("PROTON", "all")
```

```
STOREPAR("PULPROG", "zg30")  
STOREPAR("PL 1", 10.0)  
GETINT("Please enter the number of scans:", i1)  
STOREPAR("NS", i1)  
ZG  
EFP  
QUIT
```

SEE ALSO

STOREPARS - store a status parameter
STOREPARN - store a parameter to specified direction
FETCHPAR - get an acquisition, processing or output parameter

STOREPARN

NAME

STOREPAR - store a parameter to the specified direction

SYNTAX

STOREPARN(dir, par, val)

DESCRIPTION

TOPSPIN 2.1 and newer offers the macro STOREPARN. It works like STOREPAR except that it can be used for any direction of an n-dimensional dataset. STOREPARN takes three arguments:

1. the direction of the dataset
1. the name of the parameter
2. the value to be stored which can specified in two different forms:
 - as a constant
 - as the name of an AU variable

STOREPARN works on nD datasets of any dimension.

TOPSPIN 2.0 and older only supported AU parameter storage up to 3D, using the macros STOREPAR, STOREPAR1 and STOREPAR3. In TOPSPIN 2.1 and newer, these macros can still be used or they can be replaced by STOREPARN. Note that the direction specification for STOREPARN is different from STOREPAR/1/3.

For a 2D dataset:

F2 direction (acquisition direction):

STOREPAR(par, val) or STOREPARN(2, par, val)

F1 direction:

STOREPAR1(par, val) or STOREPARN(1, par, val)

For a 3D dataset:

F3 direction (acquisition direction):

STOREPAR(par, val) or STOREPARN(3, par, val)

F2 direction:

STOREPAR1(par, val) or STOREPARN(2, par, val)

F1 direction:

STOREPAR3(par, val) or STOREPARN(1, par, val)

SEE ALSO

STOREPAR - store a parameter in acquisition direction

STOREPARNS - store a status parameter to specified direction

FETCHPARN - fetch a parameter from specified direction

STOREPARS

NAME

STOREPARS - store a status parameter (acquisition and processing)

SYNTAX

STOREPARS(par, val)

DESCRIPTION

The macro STOREPARS stores the value of an AU variable into a status parameter. This AU variable can then be used in subsequent AU statements. Status parameters are set by an acquisition or processing command and describe the status of the dataset after this acquisition or processing command. If the data are now manipulated by AU statements which do not update the status parameters, these must be set explicitly with STOREPARS. For example, if you add two fid's with *addfid*, the total number of scans of the resulting dataset is not automatically updated. This must be done explicitly with STOREPARS.

The handling of the macros STOREPAR1S and STOREPAR3S is equivalent to the handling of STOREPARS.

EXAMPLE

The following AU program reads the foreground dataset, adds the fid of the next experiment number and the experiment number after that and stores the result in the foreground dataset. The number of scans of the original fid's are added and stored in the status parameter NS of the resulting dataset.

```
int expno_save;
DATASET2(name, expno+1, procno, disk, user)
DATASET3(name, expno+2, procno, disk, user)
expno_save=expno;
IEXPNO
FETCHPARS("NS", &i1)
IEXPNO
```

```
FETCHPARS("NS", &i2)
REXPNO(expno_save)
ADDFID
STOREPARS("NS", i1+i2)
QUIT
```

SEE ALSO

FETCHPARS - get a status parameter (acquisition and processing)
STOREPARNS - store a status parameter to specified direction
STOREPAR - store an acquisition, processing or output parameter

RPAR

NAME

RPAR - read a parameter set to the current AU dataset

SYNTAX

RPAR(char *parset, char *typ)

DESCRIPTION

The macro RPAR reads a parameter set to the current AU dataset. This can be a standard Bruker parameter set or a user defined parameter set which was stored with WPAR. RPAR takes two arguments:

1. the name of the parameter set
2. the type of parameters which are to be read

The second argument can be:

- *acqu* for acquisition parameters (**eda**)
- *proc* for processing parameters (**edp**)
- *outd* for output parameters (**edp**)
- *all* for acquisition, processing, plot and output parameters

EXAMPLE

The following AU program reads the standard Bruker parameter set PROTON, sets the number of scans to 1024 and runs an acquisition:

```
RPAR("PROTON", "all")  
STOREPAR("NS", 1024)  
ZG  
QUIT
```

SEE ALSO

WPAR - write the current dataset parameters to a parameter set

WPAR

NAME

WPAR - write the current dataset parameters to a parameter set

SYNTAX

WPAR(char *parset, char *typ)

DESCRIPTION

The macro WPAR writes the parameters of the current AU dataset to a parameter set. You can only write to user defined parameter sets. Bruker standard parameters sets cannot be overwritten. WPAR is typically used in AU programs to store a temporary parameter set. It takes two arguments:

1. the name of the parameter set
2. the type of parameters which are to be stored

The second argument can be:

- *acqu* for acquisition parameters (**eda**)
- *proc* for processing parameters (**edp**)
- *outd* for output parameters
- *all* for acquisition, processing, plot and output parameters

EXAMPLE

The following AU program reads the acquisition parameters of the Bruker standard parameter set PROTON, sets the number of scans, the frequency offset and time domain data size and writes the acquisition parameters to a temporary parameter set. It then performs 8 successive acquisitions with these parameters.

```
RPAR("PROTON", "all")
STOREPAR("NS", 16)
STOREPAR("O1", 766.23)
STOREPAR("TD", 8192)
WPAR("partemp", "acqu")
TIMES(8)
```

```
ZG
IEXPNO
RPAR("partemp", "acqu")
END
QUIT
```

SEE ALSO

RPAR - read a parameter set to the current AU dataset

Chapter 10

Macros for Plot Editor/autoplot

This chapter contains a description of AU macros which can be used to plot data using Plot Editor portfolios and layouts. These include macros for the creation and definition of portfolios and for plotting to the printer, to a post-script file or enhanced metafile.

TOPSPIN 2.0 and newer also offer macros for automatic creation of Plot Editor layouts. Examples are LAYOUT_OBJ_1D and LAYOUT_ADD. These are described in a separate manual (see *Help* ' *Manuals* ' **[Programming Manuals]** *Plot Layout Programming*).

AUTO PLOT

NAME

AUTO PLOT - plot the current dataset according a Plot Editor layout

SYNTAX

AUTO PLOT

DESCRIPTION

The macro AUTO PLOT plots the current dataset according to the Plot Editor layout that is defined by the parameter LAYOUT.

The Plot Editor layout can be:

- a standard layout that was delivered with TOPSPIN
- a user defined layout that was setup and stored from Plot Editor

Processing AU programs that contain the AUTO PLOT macro can be used with one of the options **a**, **e**, **h** or **t**. They cause AUTO PLOT to store the plot as a postscript file. For example, the AU program **proc_1d** can be enter as:

proc_1d - prints to the printer defined in the layout

proc_1d a - prints to a PDF file in the dataset procno

proc_1d e - also prints to a postscript file in the dataset procno

proc_1d h - also prints to a postscript file in the users home directory

proc_1d t - also prints to a postscript file in the TEMP directory

Furthermore, you can use multiple arguments, e.g.:

proc_1d a e - prints to a PDF file and to a postscript file in the dataset procno

(see also the header of the AU program **plot_to_file**)

EXAMPLE

This AU program processes the current 1D dataset and plots it according to the Plot Editor layout specified in **edp**:

EF

APK
SREF
ABS
AUTOPLOT
QUIT

SEE ALSO

AUTO_PLOT_TO_FILE, AUTO_PLOT_WITH_PORTFOLIO,
AUTO_PLOT_WITH_PORTFOLIO_TO_FILE

AUTO PLOT_TO_FILE

NAME

AUTO PLOT_TO_FILE - as AUTO PLOT but store the output into a file

SYNTAX

AUTO PLOT_TO_FILE(filename)

DESCRIPTION

The macro AUTO PLOT_TO_FILE plots the current dataset according to the Plot Editor layout defined by the parameter LAYOUT. The output is not sent to the printer but stored in the file that is specified as an argument. The argument is normally a full pathname. If it is not, the file is stored in the TOPSPIN home directory.

If the filename has the extension `.ps`, the output is stored as a postscript file. If (under Windows) it has the extension `.emf`, as in the example below, the output will be stored as an enhanced metafile.

AUTO PLOT_TO_FILE is actually a composite macro that consists of several PORTFOLIO*/AUTO PLOT* macros. This, however, is transparent to the user.

EXAMPLE

This AU program processes the current 1D dataset and plots it according to the Plot Editor layout specified in **edp**. The result is stored in an enhanced metafile.

```
EF
APK
SREF
ABS
AUTO PLOT_TO_FILE("C:/mydata.emf")
QUIT
```

SEE ALSO

AUTO PLOT, AUTO PLOT_WITH_PORTFOLIO

CREATE_PORTFOLIO

NAME

CREATE_PORTFOLIO - create a Plot Editor portfolio

SYNTAX

CREATE_PORTFOLIO(name)

DESCRIPTION

The macro CREATE_PORTFOLIO creates the Plot Editor portfolio that is specified as an argument. It takes one argument; the filename of the portfolio.

The argument is normally specified as a full pathname. If it is not, the portfolio is stored under the TOPSPIN home directory. If the specified file already exists, it is overwritten. Note that CREATE_PORTFOLIO creates the portfolio but does not insert any dataset specifications.

EXAMPLE

This AU program plots the current dataset according to the Plot Editor layout specified in *edp*. It is just a simple demonstration of the use of PORTFOLIO macros.

```
CREATE_PORTFOLIO("C:/temp/myPortfolio.por")
ADD_CURDAT_TO_PORTFOLIO
CLOSE_PORTFOLIO

AUTOPLOT_WITH_PORTFOLIO

QUIT
```

Note that this AU program does the same as the command *autoplot*.

SEE ALSO

ADD_TO_PORTFOLIO, CLOSE_PORTFOLIO

ADD_CURDAT_TO_PORTFOLIO

NAME

ADD_CURDAT_TO_PORTFOLIO - add the current dataset to the portfolio

SYNTAX

```
ADD_CURDAT_TO_PORTFOLIO
```

DESCRIPTION

The macro ADD_CURDAT_TO_PORTFOLIO adds the current dataset to the Plot Editor portfolio that was previously create with CREATE_PORTFOLIO.

EXAMPLE

This AU program plots two datasets, the current and next processing number of the current data name, according to the Plot Editor layout

```
CREATE_PORTFOLIO("C:/temp/myPortfolio.por")
ADD_CURDAT_TO_PORTFOLIO
IPROCNO
ADD_CURDAT_TO_PORTFOLIO
CLOSE_PORTFOLIO

AUTOPLOT_WITH_PORTFOLIO
QUIT
```

SEE ALSO

CREATE_PORTFOLIO, CLOSEPORTFOLIO

ADD_TO_PORTFOLIO

NAME

ADD_TO_PORTFOLIO - add the specified dataset to the portfolio

SYNTAX

ADD_TO_PORTFOLIO(disk,user, name, expno, procno)

DESCRIPTION

The macro ADD_TO_PORTFOLIO adds a dataset to the portfolio that was previously created with CREATE_PORTFOLIO. The dataset to be added is completely specified by the five arguments of ADD_TO_PORTFOLIO. Note that these arguments can be constants (values) or variables.

EXAMPLE

This AU program plot two datasets according to the TOPSPIN layout. Note that the first dataset to be plotted is the so called second dataset (**edc2**), specified by the predefined dedicated variables disk2, user2 etc.

```
CREATE_PORTFOLIO("/temp/myPortfolio.por")
GETCURDATA2
ADD_TO_PORTFOLIO(disk2, user2, name2, expno2, procno2)
ADD_TO_PORTFOLIO("C:/ts", "guest", "mydata", 1, 3)
CLOSE_PORTFOLIO

AUTOPLOT_WITH_PORTFOLIO

QUIT
```

SEE ALSO

ADD_CURDAT_TO_PORTFOLIO

CLOSE_PORTFOLIO

NAME

CLOSE_PORTFOLIO - closes the portfolio definition

SYNTAX

CLOSE_PORTFOLIO

DESCRIPTION

The macro CLOSE_PORTFOLIO closes the definition of the portfolio that was previously created with CREATE_PORTFOLIO. It must be used after the last ADD_CURDAT_TO_PORTFOLIO or ADD_TO_PORTFOLIO macro and before the first AUTO PLOT* macro.

EXAMPLE

This AU program plots the current dataset according to the TOPSPIN layout. It is just a simple demonstration of the use of PORTFOLIO macros.

```
CREATE_PORTFOLIO("C:/temp/myPortfolio.por")
ADD_CURDAT_TO_PORTFOLIO
CLOSE_PORTFOLIO

AUTO PLOT_WITH_PORTFOLIO

QUIT
```

Note that this AU program does the same as the command *autoplot*.

SEE ALSO

CREATE_PORTFOLIO, ADD_TO_PORTFOLIO

AUTO PLOT_WITH_PORTFOLIO

NAME

AUTO PLOT_WITH_PORTFOLIO - plot the dataset(s) of the current portfolio

SYNTAX

AUTO PLOT_WITH_PORTFOLIO

DESCRIPTION

The macro AUTO PLOT_WITH_PORTFOLIO plots the dataset(s) defined in the portfolio created with CREATE_PORTFOLIO according to the Plot Editor layout defined by the **edp** parameter LAYOUT.

EXAMPLE

This AU program plots the current dataset according to the TOPSPIN layout. It is just a simple demonstration of the use of PORTFOLIO macros.

```
CREATE_PORTFOLIO("C:/temp/myPortfolio.por")
ADD_CURDAT_TO_PORTFOLIO
CLOSE_PORTFOLIO

AUTO PLOT_WITH_PORTFOLIO

QUIT
```

Note that this AU program does the same as the command **autoplot**.

SEE ALSO

AUTO PLOT, AUTO PLOT_WITH_PORTFOLIO_TO_FILE

AUTO PLOT_WITH_PORTFOLIO_TO_FILE

NAME

AUTO PLOT_WITH_PORTFOLIO_TO_FILE - plot the dataset(s) of the current portfolio and store the output into a file

SYNTAX

AUTO PLOT_WITH_PORTFOLIO_TO_FILE(filename)

DESCRIPTION

The macro AUTO PLOT_WITH_PORTFOLIO_TO_FILE plots the dataset(s) defined in the Plot Editor portfolio that was previously created with CREATE_PORTFOLIO. The plot is made according to the layout defined by the parameter LAYOUT. The output is stored in the file that is specified as an argument to the macro. The argument is normally a full pathname. If it is not, the portfolio is stored under the TOPSPIN home directory.

If the filename has the extension `.ps`, as in the example below, the output will be stored as a postscript file. If (under Windows) it has the extension `.emf`, the output is stored as an enhanced metafile.

EXAMPLE

This AU program plots the current dataset according to the Plot Editor layout specified in **edp** and stores the result into a postscript file.

```
CREATE_PORTFOLIO("C:/temp/myPortfolio.por")
ADD_CURDAT_TO_PORTFOLIO
CLOSE_PORTFOLIO

AUTO PLOT_WITH_PORTFOLIO_TO_FILE("/users/guest/mydata.ps")

QUIT
```

SEE ALSO AUTO

PLOT_WITH_PORTFOLIO, AUTO PLOT_TO_FILE

Chapter 11

Macros prompting the user for input.

This chapter contains a description of all AU macros which can be used to prompt the user for input and store the input into an AU variable. Different macros are available for requesting integer, float, double or text values.

GETINT

NAME

GETINT - prompt the user to enter an integer value

SYNTAX

```
GETINT("Please enter an integer value : ", i1)
```

DESCRIPTION

The macro GETINT prompts the user to enter an integer value and stores this value into an integer variable. It can be used for various purposes, for example to set the value of a TOPSPIN (integer) parameter or to specify the number of cycles in an AU program loop. GETINT takes two arguments:

1. a text string which prompts the user to enter an integer value
2. an integer variable into which the value is stored

EXAMPLE

The following AU program prompts the user for the number of scans per acquisition and the number of experiments to be done:

```
GETINT("Please enter the number of scans:", i1)
GETINT("Please enter the number of experiments:", i2)
TIMES(i2)
  STOREPAR("NS", i1)
  ZG
  IEXPNO
END
QUIT
```

SEE ALSO

GETFLOAT - prompt the user to enter a float value
GETDOUBLE - prompt the user to enter a double value
GETSTRING - prompt the user to enter a text string

GETFLOAT, GETDOUBLE

NAME

GETFLOAT - prompt the user to enter a float value

GETDOUBLE - prompt the user to enter a double value

SYNTAX

GETFLOAT(text, f1)

GETDOUBLE(text, d1)

DESCRIPTION

The macro GETFLOAT prompts the user to enter a float value and stores this value into a float AU variable. It is used to get the value for a TOPSPIN (float) parameter which can then be stored with STOREPAR. GETFLOAT takes 2 arguments:

1. a text string which prompts the user to enter an float value
2. the float variable into which the value is store

The description for GETDOUBLE is equivalent, except that it is used for TOPSPIN (double) parameters.

EXAMPLE

The following AU program prompts the user for the *frequency offset* and *Gaussian broadening*, stores these values into the parameters O1 and GB respectively and then runs an acquisition, Gaussian multiplication and Fourier transform:

```
GETDOUBLE("Please enter the frequency offset:", d1)
STOREPAR("o1", d1);
GETFLOAT("Please enter the Gaussian broadening:", f1)
STOREPAR("GB", f1)
ZG
GM
FT
QUIT
```

SEE ALSO

GETINT - prompt the user to enter an integer value

GETSTRING - prompt the user to enter a text string

GETSTRING

NAME

GETSTRING - prompt the user to enter a text string

SYNTAX

GETSTRING(text, cval)

DESCRIPTION

The macro GETSTRING prompts the user to enter a text string which is then stored into an AU variable. It can be used for various purposes, for example to ask the user a question or prompt the user to enter a name. GETINT takes two arguments:

1. a text string which prompts the user to enter a text string
2. the character-string variable into which the value is stored

EXAMPLE

The following AU program asks the user if an integration must be done and, if yes, which intrng file must be used. Then the integrals are calculated and listed:

```
char answer[8];  
  
GETSTRING("Do you want to do an integration (yes/no)?", answer)  
if ( !strcmp(answer,"yes") )  
{  
    GETSTRING("Which intrng file must be used?", text)  
    RMISC("intrng", text)  
    LI  
}  
QUIT
```

SEE ALSO

GETINT - prompt the user to enter an integer value
GETFLOAT - prompt the user to enter a float value
GETDOUBLE - prompt the user to enter a double value

Chapter 12

Bruker library functions

This chapter contains a description of various C functions which are available as part of Bruker libraries. Several of them concern the handling of dataset lists or directory lists. You can, for instance, get a list of filenames, display it, select a file from the list and then copy the file to a different directory. The functions described in this chapter are particularly useful for files located in the directories `<tshome>/conf` and `<tshome>/exp`. For copying datasets, we recommend to use the macros described in Chapter 5.

CalcExpTime, PrintExpTime

NAME

`CalcExpTime` - calculate the experiment time for the current experiment
`PrintExpTime` - print the experiment time for the current experiment

SYNTAX

```
static void PrintExpTime();  
  
int CalcExpTime ();  
void PrintExpTime (int exptime, int expno);  
  
#include<inc/exptutil>
```

DESCRIPTION

The function `CalcExpTime` calculates the experiment time for the current experiment. The return value is the experiment time in seconds. The function `PrintExpTime` can be used to print the experiment time in the form "days hours minutes seconds".

EXAMPLE

The following AU program calculates and prints the experiment time of a sequence of 10 experiments starting with the foreground dataset.

```
static void PrintExpTime();  
  
TIMES(10)  
  PrintExpTime (CalcExpTime (), loopcount1);  
  IEXPNO  
END  
QUIT  
#include<inc/exptutil>
```

Note that the declaration of `PrintExpTime` must appear at the beginning of the AU program and the `#include` statement at the end of the AU program.

SEE ALSO

multiexpt - a standard Bruker library AU program

ChecksumFile

NAME

`ChecksumFile` - creates a checksum of a data file

SYNTAX

`ChecksumFile(filnam, 0, dest, 0, verb, bytord, dtyp, dim, siz0, siz, xdim)`

DESCRIPTION

The function `ChecksumFile` generates a checksum of a data file. The output consist of a checksum preceded by hash type and data sizes, e.g:

```
data hash MD5: 512 * 256
```

```
16 A5 E9 14 FB 66 8B 48 09 8B E3 CA 86 D2 68 A2
```

which stored in a destination character string. The input data file can be a TOPSPIN raw or processed data file or any other integer data file. The data size, storage mode and dimensionality must be specified as arguments.

The arguments of the function have the following meaning:

- `const char* filnam`
Full pathname of the input data file.
- `char* dest`
Destination character string for function output. Must at least be 128 byte.
- `int verb`
Verbose error if the input file does not exist (0=no, 1=yes)
- `int bytord`
Byteorder of the input data (parameter `BYTORDA` for raw data or `BYTORDP` for processed data).
- `int dtyp`
Data type of the input data (parameter `DTYPA` for raw data and `DTYPP` for processed data)
- `int dim`
Data dimensionality (1 for 1D data, 2 for 2D data, ... etc.)

- `int siz0`
For raw data, `siz0` must be set to TD in the acquisition direction. For processed data, `siz0` must be set to SI in the first direction.
 - `const int* siz`
Array of sizes. For processed data, `siz` must be set to SI in the successive directions. For example, for 2D data, `siz = (F2-SI, F1-SI)`. For raw data, `siz` must be set to TD in the successive directions.
- Attention: in the acquisition direction, `siz` must be TD rounded to the next multiple of 256.
- `const int* xdim`
For raw data, `xdim = siz`. For processed data, `xdim` is the array of submatrix sizes of successive directions. For example, for 2D processed data `xdim = (F2-XDIM, F1-XDIM)`.

ChecksumFile can have the following return values:

- > 0 : successful operation
- = 0 : parameter inconsistency or I/O problems
- < 0 : all other cases

The return value can be used as an argument of the function CheckSumError() which generates an error string.

The output of CheckSumFile can be appended to the audit file with the function AuditAppend as shown in the example below.

EXAMPLE

The following AU statements will generate a data checksum of the current processed 2D data and store it in the current data auditp.txt file. It is part of the Bruker AU program *xfshear*.

```
int bytordp, dtyp, size[2], xdim[2];
char name2r[PATH_MAX], nameaudit[PATH_MAX], audittext[512];
char* auditp = audittext + sprintf(audittext, "processing description");
FETCHPARS("BYTORDP",&bytordp)
FETCHPARS("DTYPP",&dtyp)
FETCHPARS("SI",&size[0])
```

```
FETCHPARS("XDIM", &xdim[0])
FETCHPAR1S("SI",&size[1])
FETCHPAR1S("XDIM", &xdim[1])

sprintf(name2rr,"%s/data/%s/nmr/%s/%d/pdata/%d/2rr", disk, user,
name, expno, procno);

if (CheckSumFile(name2rr, 0, auditp, 0, 1, bytordp, dtypp, 2, size[0],
size, xdim) > 0)
{
    sprintf(nameaudit,"%s/data/%s/nmr/%s/%d/pdata/%d/au-
ditp.txt",disk, user, name, expno, procno);
    AuditAppend(nameaudit, audittext);
}

QUIT
```

SEE ALSO

AuditAppend, AuditCreate

AuditAppend

NAME

`AuditAppend` - appends a new entry to an existing audit file

SYNTAX

```
AuditAppend(const char* auditn, const char* what)
```

DESCRIPTION

The function `AuditAppend` appends a new entry to an existing audit file. An audit file entry consists of the following fields:

NUMBER, WHEN, WHO, WHERE, PROCESS, VERSION, WHAT

All of these are automatically set by `AuditAppend`, except for the `WHAT` field which is specified as the second argument. It can be any character string.

Note that `AuditAppend` does not create an audit file if this does not exist yet.

EXAMPLE

See the example of the function `ChecksumFile`.

SEE ALSO

`ChecksumFile`, `AuditCreate`

AuditCreate

NAME

`AuditCreate` - create a new audit file

SYNTAX

```
AuditCreate(const char* auditn, const char* what)
```

DESCRIPTION

The function `AuditCreate` creates a new audit file with a single entry. This is, for example, useful when new data are created. An audit trail entry consists of the following fields:

NUMBER, WHEN, WHO, WHERE, PROCESS, VERSION, WHAT

All of these are automatically set, except for the WHAT field which is specified as the second argument.

Caution: `AuditCreate` overwrites a possibly existing audit file.

EXAMPLE

Please look at the Bruker AU program ***split3d*** for an example of using `AuditCreate`.

SEE ALSO

`ChecksumFile`, `AuditAppend`

FileSelect

NAME

`FileSelect` - display a list of directory entries and allow to select entries

SYNTAX

```
FileSelect(sourcedir, flist, &num, type);
```

DESCRIPTION

The function `FileSelect` opens a directory, shows a list of all file and directory entries and allows you to select one or more entries. The list is stored in a 2 dimensional character-string variable which can be evaluated by subsequent AU statements. `FileSelect` takes four arguments:

1. the source directory
2. the variable into which the list is stored
3. the variable into which the number of selected entries is stored
4. a flag which determines whether files (0) or directories (1) are listed.

`FileSelect` replaces the functions `uxselect` and `getdir` which were used by TOPSPIN's predecessor XWIN-NMR.

EXAMPLES

The following AU program will make a list of all shim files and will display this list in a selection window. If an entry is selected, then the corresponding shim file is read with the macro RSH. If no entries were found or selected, the AU program aborts.

```
char sourcedir[200];
int num= 100;
char flist[128][128];

sprintf(sourcedir, "%s/nmr/lists/bsms", PathXWinNMRExpStan());
if (i1 = FileSelect(sourcedir, flist, &num, 0) < 0)
{
    Proc_err(DEF_ERR_OPT, "Error: No shim files selected!");
}
```

```
    }  
    else  
    {  
        RSH(flist[0])  
    }  
    QUIT
```

getdir

NAME

`getdir` - get all file names and/or directory names within a directory

SYNTAX

```
int getdir (char *directory, char ***filelist, const char *match-code);
```

DESCRIPTION

The function `getdir` opens a directory and gets all file and directory names in that directory. This list is stored in a 2 dimensional character-string variable which can be evaluated by subsequent AU statements. The list can be limited by specifying a match-code; only names matching this string are entered into the list. `getdir` takes three arguments:

1. the source directory
2. the variable into which the list of names is stored
3. the match-code; an arbitrary string of characters

The third argument can also be `"/files"` to get all files but not directories, or `"/dir"` to get all directories but not files.

The return value of `getdir` is the number of successfully matched file names and/or directory names.

`getdir` internally allocates memory for the list of names. Officially, you must free this memory with the Bruker library function `freedir`. In practice, however, you can omit `freedir` because all memory allocated by the AU program is automatically freed when the AU program finishes.

EXAMPLES

The following AU statements will create a list of experiment directories from an TOPSPIN dataset. All entries are returned because no match-code was specified.

```
char sourcedir[200], **listfile;  
sprintf (sourcedir, "%s/data/%s/%s/%s/", disk, user, type, name);  
i1 = getdir (sourcedir, &listfile, NULL);
```

The following AU statements will create a list of shim files starting with the letters a to p from the bsms directory.

```
char sourcedir[200], **listfile;
sprintf(sourcedir, "%s/nmr/lists/bsms", PathXWinNMRExpStan());
i1 = getdir (sourcedir,&listfile,"[a-p]*");
```

The following AU statement will create a list of all files but not directories from the users home directory.

```
i1 = getdir (PathSystemHome(),&listfile,"/files");
```

The following AU statement will return a list of all directory names from the users home directory.

```
i1 = getdir (PathSystemHome(),&listfile,"/dir");
```

SEE ALSO

`freedir` - free memory allocated by `getdir`

freedir

NAME

`freedir` - free memory allocated by `getdir`

SYNTAX

```
void freedir (char **listfile);
```

DESCRIPTION

The function `freedir` frees the memory that was allocated by a `getdir` function call.

EXAMPLE

See the example under the function `FileSelect`.

SEE ALSO

`getdir` - get all file names and/or directory names within a directory
`fileselect` - display a list from which an entry can be selected by mouse-click

dircp, dircp_err

NAME

`dircp` - copy a file

`dircp_err` - return the error message that corresponds to the error return value of a `dircp` function call

SYNTAX

```
dircp (char *sourcefile, char *targetfile);  
char *dircp_err (int return-value);
```

DESCRIPTION

The function `dircp` copies the sourcefile into the targetfile. If the targetfile exists, it will be overwritten. A negative number is returned if copying was not possible. The function `dircp_err` will return the corresponding C error message. A return value of 0 indicates successful execution.

EXAMPLE

The following AU program copies the `title` file of the foreground dataset to the users home directory.

```
char sourcefile[200], targetfile[200];  
  
sprintf (targetfile, "%s/title", PathSystemHome());  
if ( (i1 = dircp (PROCPATH("title"),targetfile)) < 0 )  
    Proc_err (DEF_ERR_OPT, dircp_err(i1));  
QUIT
```

Note that `PROCPATH` uses a static buffer for building the pathname, which means it cannot be used to built more then one pathname at a time, e.g. it cannot be used in both `dircp` arguments.

gethighest

NAME

`gethighest` - return the next highest unused experiment number of a dataset

SYNTAX

```
int gethighest (char *directory);  
#include <inc/sysutil>
```

DESCRIPTION

The function `gethighest` scans a directory for all subdirectories whose name is a number and returns then returns the next highest unused number. `gethighest` is typically used to scan a dataset name directory of a TOPSPIN dataset. In that case, it returns the highest unused experiment number. If, for example, the highest used experiment number is 56, the function will return the value 57. The function can also be used to return the highest unused processing number of a dataset.

EXAMPLE

The following AU program will copy the current TOPSPIN experiment into the next highest unused experiment dataset.

```
(void) sprintf (text,"%s/data/%s/nmr/%s",disk,user,name);  
i1 = gethighest (text);  
WRA(i1)  
QUIT  
#include <inc/sysutil>
```

Note that the `#include` statement must be included at the end of the AU program.

getParfileDirForRead

NAME

getParfileDirForRead - determine the pathname of a list file to be read

SYNTAX

```
int getParfileDirForRead (const char *filename, const char *key, char
*path);
```

DESCRIPTION

The function `getParfileDirForRead` determines the pathname of a list file to be read. The function has three arguments:

- the filename of the list file
- the type (key) of the list file: `PP_DIRS`, `VD_DIRS` etc. (see Figure 12.1)
- the pathname of the list file

The third argument contains the result of the function. For determining this pathname, the function searches for the specified filename in all source directories that are set up for the specified list type, for the current user. The first source directory in which the file is found determines the output pathname. To view or change the list of source directories:

1. Click *Options* ' *Preferences* [*set*]
2. Click *Miscellaneous* in the left part of the dialog box.
3. Click the *Change* button of the *entry Manage source directories...*

The functions `getParfileDirForRead` and `getParfileDirForWrite` are only implemented in TOPSPIN 2.1 and newer. The replace the functions `getstan` and `PathXWinNMR*`.

EXAMPLE

The following AU statements are an example of the usage of the function `getParfileDirForRead`. They are part of the AU program *proc_intrng*.

```
char intrngfilePathname[PATH_MAX];
if (getParfileDirForRead("testrng", INTRNG_DIRS, intrngfilePathname)
```

```
< 0)
{
Proc_err(DEF_ERR_OPT, "testrng: %s", intrngfilePathname);
ABORT
}
RMISC("intrng", intrngfilePathname)
```

SEE ALSO

getParfileDirForWrite - Determines pathname of list file to be written.

PP_DIRS	pulse programs
CPD_DIRS	cpd programs
MAC_DIRS	macros
PY_DIRS	python programs
GP_DIRS	gradient programs
SHAPE_DIRS	shape files
SP_DIRS	shape lists
AU_DIRS	AU programs
PAR_DIRS	parameter sets
VD_DIRS	delays
VP_DIRS	pulses
VC_DIRS	loop counters
VT_DIRS	temperatures
VA_DIRS	amplitudes
F1_DIRS	frequencychannel 1
DS_DIRS	data sets
SCL_DIRS	solvent scaling region
PHASE_DIRS	phase programs
INTRNG_DIRS	intrng files
PEAKRNG_DIRS	peakrng files
BASLPNTS_DIRS	baslpnts files
BASE_INFO_DIRS	base_info files
PEAKLIST_DIRS	peaklist files
CLEVELS_DIRS	clevels files
REG_DIRS	region files
INT2DRNG_DIRS	int2drng files

Figure 12.1

getParfileDirForWrite

NAME

getParfileDirForWrite - determine the pathname of a list file to be written

SYNTAX

```
int getParfileDirForWrite (const char *filename, const char *key, char
*path);
```

DESCRIPTION

The function `getParfileDirForWrite` determines the pathname of a list file to be read. The function has three arguments:

- the filename of the list file
- the type (key) of the list file: `PP_DIRS`, `VD_DIRS` etc. (see Figure 12.1 at the description of `getParfileDirForRead`)
- the pathname of the list file

The third argument contains the result of the function. For determining this pathname, the function searches through all source directories that are set up for the specified list type, for the current user. The first source directory that actually exists on disk ¹, determines the output pathname. To view or change the list of source directories:

1. Click *Options* ' *Preferences* [*set*]
2. Click *Miscellaneous* in the left part of the dialog box.
3. Click the *Change* button of the *entry Manage source directories...*

The functions `getParfileDirForRead` and `getParfileDirForWrite` are only implemented in TOPSPIN 2.1 and newer. The replace the functions `getstan` and `PathXWinNMR*`.

EXAMPLE

The following AU statements are an example of the usage of the function

1. Usually the first specified source directory

getParfileDirForWrite. They are part of the AU program **sysgenpar**.

```
char vspath[PATH_MAX];
```

```
if (getParfileDirForWrite("syst1list", VD_DIRS, vspath) < 0)
```

```
{
```

```
    Proc_err(DEF_ERR_OPT, "syst1list: %s", vspath);
```

```
    ABORT
```

```
}
```

SEE ALSO

getParfileDirForRead - Determines pathname of list file to be read.

getstan

NAME

`getstan` - return the pathname to the user's current experiment directory

SYNTAX

```
char *getstan (char *pathname, const char *subdir);
```

DESCRIPTION

The function `getstan` returns the pathname to the user's current experiment directory. The returned pathname can be concatenated with a known subdirectory pathname as a part of the same `getstan` function call.

Please note: in TOPSPIN 2.1 and newer the functions `getstan` and `PathXWinNMR*` have become obsolete and can be replaced by the functions `getParfileDirForRead` and `getParfileDirForWrite`.

EXAMPLE

The following AU program will get the pulse program of the current AU dataset. It will then prompt the user to confirm the name of the pulse program or to enter a new name. Finally, the pulse program will be shown in a TOPSPIN window.

```
char pulprog[80];  
  
FETCHPAR("PULPROG",pulprog)  
GETSTRING ("Enter the name of the pulse program: ",pulprog);  
(void) sprintf (text,"%s/%s",getstan (NULL,"lists/pp"),pulprog);  
showfile (text);  
QUIT
```

NOTE

In the above example, the function call `getstan (NULL,"lists/pp")` returns the pathname `/<tshome>/exp/stan/nmr/lists/pp`. The function call `getstan(NULL,NULL)` returns `/<tshome>/exp/stan/nmr/`.

SEE ALSO

`PathXWinNMR*` - a class of functions which return pathnames to certain TOPSPIN directories

`getParfileDirForRead` - Determines pathname of list file to be read.

`getParfileDirForWrite` - Determines pathname of list file to be written.

getxwinvers

NAME

`getxwinvers` - return the current version and patchlevel of TOPSPIN

SYNTAX

```
int getxwinvers (char *curversion);  
#include <inc/sysutil>
```

DESCRIPTION

The function `getxwinvers` returns the version and patchlevel of the currently running TOPSPIN program into the variable `curversion`. This variable can then be printed out.

EXAMPLE

The following AU program prints the current version and patchlevel in the status line of TOPSPIN.

```
char curversion[80];  
i1 = getxwinvers(curversion);  
show_status (curversion);  
QUIT  
#include <inc/sysutil>
```

Note that the `#include` statement must be included at the end of the AU program.

mkudir

NAME

`mkudir` - create a complete directory path

SYNTAX

```
int mkudir (char *directory);
```

DESCRIPTION

The function `mkudir` scans the specified directory for the last `/`. Then it checks recursively for the existence of all components of the directory path and creates them if necessary. The function returns `-1` if an error occurs, otherwise `0`.

If the full pathname is to be created, then the directory must end with a `/` (see the example below). Possible characters behind the last slash are discarded.

EXAMPLE

The following AU program will create a dataset directory tree which has an experiment number one higher than the current foreground dataset.

```
(void) sprintf (text, "%s/data/%s/nmr/%s/%d/pdata/%d",
                disk, user, name, expno+1, procno);
if (mkudir(text) < 0)
    Proc_err (DEF_ERR_OPT, "could not create :\n%s", text);
QUIT
```

PathXWinNMR

NAME

PathXWinNMR - a class of functions which return pathnames to certain TOPSPIN directories

SYNTAX

```
char *PathXWinNMRConf ();
char *PathXWinNMRCurDir ();
char *PathXWinNMRDotXWinNMR ();
char *PathXWinNMRExp ();
char *PathXWinNMRPlot ();
char *PathXWinNMRProg ();
```

DESCRIPTION

The above functions return pathnames to certain TOPSPIN mostly subdirectories of the TOPSPIN directory <tshome>. For a standard installation, <tshome> is:

- on LINUX systems: /opt/topspin
- on Windows systems: C:\Bruker

Please note: in TOPSPIN 2.1 and newer the functions `getstan` and `PathXWinNMR*` have become obsolete and can be replaced by the functions `getParfileDirForRead` and `getParfileDirForWrite`.

For a user-defined installation, <tshome> can be any directory. The following table lists the directory pathnames returned by the above functions. For examples, please check the Bruker AU program library.

<code>char *PathXWinNMRConf</code>	: returns <tshome>/conf
<code>char *PathXWinNMRCurDir</code>	: returns <tshome>/prog/curdir
<code>char *PathXWinNMRDotXWinNMR</code>	: returns \$HOME/.xwinnmr-host-name
<code>char *PathXWinNMRExp</code>	: returns <tshome>/exp
<code>char *PathXWinNMRPlot</code>	: returns <tshome>/plot

`char *PathXWinNMRProg` : returns /<tshome>/prog

SEE ALSO

`getParfileDirForRead` - Determines pathname of list file to be read.

`getParfileDirForWrite` - Determines pathname of list file to be written.

pow_next

NAME

`pow_next` - round to the next larger power of two

SYNTAX

```
int pow_next (int i1);  
#include <inc/sysutil>
```

DESCRIPTION

The function `pow_next` takes `i1` and rounds it to the next larger integer value which is a power of two. The return value of the function is this power of two value. The function has no error handling. If `i1` is smaller than 1, then the function will return 1.

EXAMPLE

The following AU program will return 8192 in `i2` because this is the next larger number (compared to `i1`) which is a power of two.

```
i1 = 7000;  
i2 = pow_next(i1);  
QUIT  
#include <inc/sysutil>
```

Note that the `#include` statement must be included at the end of the AU program.

Proc_err

NAME

`Proc_err` - show a error message in a TOPSPIN dialog window

SYNTAX

```
int Proc_err (int flag, char *format);
```

```
int Proc_err (int flag, char *format, varargs);
```

DESCRIPTION

The function `Proc_err` can be used to construct a error message which will be displayed in a TOPSPIN dialog window. The function takes two or three arguments:

1. a flag which determines the type and the number (2 or 3) of buttons in the error window.
2. the error message to be displayed. If this argument contains %d, %f, or %s statements, then `Proc_err` needs a third argument which provides the corresponding variables.
3. variables who's values replace the corresponding %d, %f, or %s statements of the second argument.

The first argument (flag) can have the following values:

- `DEF_ERR_OPT`
The error window has one button (*OK*). The AU programs holds until the user clicks *OK*.
- `INFO_OPT`
The error window has one button (Seen). The AU program continues but the error window remains on the screen until it is cleared by another error window or the user clicks *Seen*.
- `QUESTION_OPT`
The error window has two buttons, OK and CANCEL. `Proc_err` returns `ERR_OK` (0) if the *OK* button is clicked and `ERR_CANCEL` (-1) if the *CANCEL* button is clicked. The return value is normally used by subsequent control statements to decide whether or not to continue the AU program.

Note that the message in the `Proc_err` window is constructed in the same way as the C function `sprintf` constructs its strings.

EXAMPLE

The following examples show several possibilities of constructing error messages for the `Proc_err` function call.

Example for `DEF_ERR_OPT` :

```
(void) sprintf (text, "%s/data/%s/nmr/%s/%d/pdata/%d",
                disk,user,name,expno+1,procno);
Proc_err (DEF_ERR_OPT, "could not create :\n%s",text);
```

Example for `QUESTION_OPT` :

```
i1 = Proc_err(QUESTION_OPT,"Continue with the AU program ?\n\
Click OK to continue, click cancel stop");
if ( i1 == ERR_OK )
{
    /* Further AU statements */
}
if ( i1 == ERR_CANCEL )
{
    ABORT
}
```

Example for `INFO_OPT` :

```
i1 = 7;
i2 = 5;
Proc_err(INFO_OPT,"%d is bigger than %d",i1,i2);
```

SEE ALSO

`Show_status` - show a string in the status line of TOPSPIN

All AU programs from the Bruker AU program library which use `Proc_err`

Show_status

NAME

Show_status - show a string in the TOPSPIN status line

SYNTAX

```
void Show_status (char *text);
```

DESCRIPTION

The function Show_status displays the specified text in the TOPSPIN status line. This function can be used as an alternative to the Proc_err function. One difference to Proc_err is that there is no window that needs to be acknowledged.

EXAMPLE

The following AU program will display the line "The AU program test has started" in the status line of TOPSPIN:

```
(void) strcpy(text,"The AU program test has started");  
Show_status (text);  
QUIT
```

SEE ALSO

Proc_err - show a message in a TOPSPIN dialog window

showfile

NAME

`showfile` - show the contents of a file in a TOPSPIN window

SYNTAX

```
int showfile (char *file);
```

DESCRIPTION

The function `showfile` reads the specified file and displays it in a TOPSPIN window. This display is a read-only display, so the file cannot be changed.

EXAMPLE

The following AU program will show the title file in a TOPSPIN window.

```
(void) sprintf (text, "%s/data/%s/nmr/%s/%d/pdata/%d/title",  
disk,user,name,expno,procno);  
i1 = showfile (text);  
QUIT
```

ssleep

NAME

`ssleep` - pause in an AU program for a certain number of seconds

SYNTAX

```
int ssleep (int seconds);
```

DESCRIPTION

The function `ssleep` will cause the AU program to wait with the execution of the next statement until the specified number of seconds has elapsed.

EXAMPLE

The following AU program will wait for 3 minutes before it resumes execution.

```
i1 = ssleep (180);  
EFP  
QUIT
```

SEE ALSO

`WAIT_UNTIL` - hold the AU program until the specified date and time

unlinkpr

NAME

`unlinkpr` - delete all processed data files (1r, 1i, 2rr, 2ii etc.) of a dataset

SYNTAX

```
int unlinkpr (char *directory);  
#include <inc/sysutil>
```

DESCRIPTION

The function `unlinkpr` deletes all processed data files (1r, 1i, 2rr, 2ii, 2ri, 2ir, dsp, dsp.hdr, dsp_low) in the specified dataset directory. There is no error check whether the files could be deleted; the return value of the function is always 0 and can be ignored.

EXAMPLE

The following AU program will delete the processed data files of the foreground dataset.

```
(void) sprintf (text, "%s/data/%s/nmr/%s/%d/pdata/%d",  
disk, user, name, expno, procno);  
i1 = unlinkpr (text);  
QUIT  
#include <inc/sysutil>
```

Note that the `#include` statement must be included at the end of the AU program.

Chapter 13

List of Bruker AU programs

13.1 Short description of all Bruker AU programs

This chapter contains a list with the names and short-descriptions of all Bruker library AU programs.

abs2.water	Performs an F2 baseline correction on a 2D dataset left and right of the water peak.
abs2D	Performs a baseline correction on a 2D dataset in both dimensions.
acqu_fid_ser	Acquires a single FID of the current 2D experiment and replaces the old fid in the ser file.
acqulist	Set up and start acquisitions using f1, f2, f3, vt, vc, vd, vp lists.
amplstab	Calculates the amplitude stability based on a peaklist file.
angle	Perform multiple acquisitions and ft's. This program is particularly interesting when you want to adjust the magic angle for MAS type experiments.
asclev	Converts the level file in the current processed data directory to ASCII and writes it to the file.
atpplot	Create a plot and a pdf file within ATP. Execute when the ATP button 'Print actual spectrum' is clicked.
au_cp	Acquisition with adjustment of decoupling power to acquisition time.
au_get1d	Acquire sweep width optimized 1D spectra.
au_getlcosy	Acquire sweep width optimized COSY spectra.
au_getlinv	Acquire sweep width optimized 2D inverse spectra.
au_getlxhco	Acquire sweep width optimized XH correlated spectra.
au_mult	AU program for C13 multiplicity analysis.
au_noediff	noe difference spectroscopy using different expnos.
au_noemult	noe difference spectroscopy with multiple irradiation points for each multiplet using different expnos.
au_water	Acquire water-suppression spectra for use in foreground (xau,xaua).
au_watersc	Acquire water-suppression spectra for use in automation, e.g., with sample changer.
au_zg	General AU program for data acquisition.

au_zg135	Acquire DEPT135 type spectra.
au_zgcosy	Acquire COSY type spectra.
au_zgglp	Automatic data evaluation according to GLP standards. This AU program takes O1, SW and O2 as arguments and then works like au_zg.
au_zgnr	Acquisition with rotation switched off.
au_zgonly	General AU program for data acquisition.
au_zgsino	Acquisition with signal to noise break up.
au_zgte	Acquisition with temperature setting.
aunmp_tojdx	Used in LIMS automation to process data. First, AUNMP is executed, then, if specified, the command given on the command line.
autoflist	Automatic generation of a frequency list for the peaks in the plot region of the spectrum.
autot1	Automatic processing of a 2D T1/T2 experiment with subsequent T1/T2 calculation.
bsms_exam	Example AU program which shows how to use low level functions to read or write BSMS parameters.
bsms_getlock	Read current LockLevel from BSMS Unit.
butselau	AU program for selective experiments in bnmr .
buttonau	AU program for basic experiments in bnmr .
butsel90	AU program for calibrating selective 90 pulse bnmr .
calcphhomo	Calculate the phase correction for the F2 and F1 dimension of homonuclear 2D experiments.
calcphinv	Calculate the phase correction for the F1 dimension in HMQC/HSQC type experiments.
calcplen	Calculate the pulse length according to the power level.
calcpowlev	Calculate the power level according to the pulse length.

calctemp	Calculate the temperature in the probe using the chemical shift difference between the aliphatic and OH protons.
calfun	Calculates an FID from an arbitrary function. This AU program is especially useful when you want to create a user defined window function for the 'uwm' command.
calibo1p1o3	Calibrate O1 and P1 for H2O samples.
check-vtu	Updates TE variable from the actual temperature and store it in a separate file (edte) in the dataset.
clev	Automatically calculate levels for 2D data.
clspec	Cleans spectra from the effects of solvent suppression. Multiple Regions can filtered or deleted from the spectrum. Entries for these regions can be deleted from peak lists and integrals. Regions are requested interactively.
coiltemp	Read the Shim Coil Temperature.
convbin2asc	Writes a 1D spectrum, with or without imaginary data points, into a file in ASCII table format. Each line in the file corresponds to one data point. The resulting file, named ascii-spec.txt, can be used to import a 1D spectrum into third party software, like Matlab.
convfidtoasc	AU program to convert an fid data file into an ascii table containing point number and intensity values. The output file is stored in the same directory as the fid. It can be used for calculations in spread-sheet programs or other third-party software.
convto1d	Converts a 2D spectrum to 1D format.
covariance	Processes data according to Covariance NMR
dcorr	Enables/disables the automatic DC-offset correction procedure of the software in case of a DRU installed.
decon_t1	Automatic deconvolution of a 2D T1/T2 experiment.
deptcyc	Creates 3 DEPT experiments from 13C experiment with CPD and then performs multiple cycles of NS scans (times 2 for DEPT90).

depthalt	Halt "deptyc" AU program.
diffe	Calculate the difference spectra between expnos.
diffp	Calculate the difference spectra between procnos.
dosy	Setup for diffusion/DOSY experiments linear gradient amplitude ramp.
elim_ints	Eliminates regions from the intrng file that contain the solvent and/or reference signals. The result will be an intrng file where the integral trails have a more reasonable scaling and smaller integrals are better resolved.
f1ref	Corrects the referencing in F1 for inverse type experiments.
fidadd	Add up FID's in incremented expno's.
fidtoser	Writes a number of fids that are stored under the same NAME and incremental EXPNOs to a ser file.
getphsum	Reads the total phase values from the status parameters and stores them back to the actual parameters.
gifadosy	Gifa starter AU program.
goalternate	Acquire alternated X/Y measurements. N averages are acquired alternatingly in two experiments.
graderror	Shows error messages generated by the gradshim gradient shimming procedure.
gontp	Starts an ntp test program.
gradratio	Calculates gradient ratios for common inverse gradient pulse programs.
gradratiogs	Calculates gradient ratios for common inverse gradient pulse programs.
gradshimau	Start gradshim gradient shimming procedure.
gsau	Program to start the gradshim gradient shimming procedure
gsel_setup	AU program to determine the transmitter offset for 1H selective gradient shimming using 1H as observe nucleus.

heater	Switch the heater on/off.
humpcal	Performs the 'hump test'. Measures the width of a peak at 0.55% and 0.11% of its signal height.
hwcal	Calculate the width of a peak at half height.
iexpno	Program to change to a new experiment number.
ift3d	Inverse Fourier transform of 3 dimensional data.
interleave	Perform interleaved acquisitions
ilhalt	Stop an interleaved acquisition which was started with the AU program interleave.
jconv_aufx	Converts Jeol FX data in a loop. The data must be stored with increasing extensions like proton.1, proton.2, ... etc.
listall_au	Scans all AU programs and extracts the name and the short description. This information is then copied into the file listall in your home directory. This list corresponds to the list you are currently reading.
loadshimZ	Reads the on-axis shim values from disk and loads them to the BSMS.
lock_off	Switch off the lock to start data acquisition on the lock channel.
lock_on	Switch on the lock if it has been disabled.
loopadj	Parameter optimization au program which calculates the lock parameters loop filter, loop gain and loop time for optimal long-time stability after adjusting lock phase and lock gain to optimal.
make2d	Create a new 2D dataset from the current 1D dataset. Can be used for 2D spectroscopy and relaxation experiments. F2 parameters are copied from the 1D data, F1 parameters are set to reasonable values.
mkflist	Automatically generates a frequency list file.
mulabel	Processing AU program for determination of ^{13}C multiplicity.

multanal	Processing AU program for determination of ^{13}C multiplicity.
multi_decon	Automatic deconvolution of a series of 1D spectra with AI calibration.
multi_integ	Automatic integration of a series of 1D spectra with AI calibration.
multi_integ2	Automatic integration of a series of 1D spectra with calibration of the integral values.
multi_integ3	Automatic integration of a series of 1D spectra with AI calibration. The output is written in a format suitable for import in excel or similar desktop publishing programs.
multi_zgvd	Performs multiple acquisitions on increasing expnos with delays that are read from a vlist file. Alternatively, a fixed delay can be entered.
multi_zgvt	Performs multiple acquisitions on increasing expnos with temperatures that are read from a vlist file.
multicmd	Performs multiple commands on increasing expnos.
multicom	Executes a TOPSPIN command in increasing expnos.
multicyc	Cycles through a series of acquisitions of increasing expnos.
multiefp	Performs multiple "efp" on increasing expnos.
multiexpt	Calculates experimental time for multizg.
multifp	Performs multiple "fp" on increasing expnos.
multiftapk	Performs multiple "ft;apk" on increasing expnos.
multihalt	Halt "multicyc" AU program.
multimas	Performs multiple MAS experiments on increasing expnos.
multiptom	Executes a TOPSPIN command in increasing procnos.
multiwinpro	Performs multiple processing on increasing expnos. The program asks for the window function and its parameters.
multixfb	Performs multiple "xfb" on increasing expnos.

multizg	Performs multiple acquisitions on increasing expnos.
noediff	noe difference spectroscopy using different expnos.
noeflist	Automatic generation of a frequency list with the peaks from the current plot region for noe.
noemult	noe difference spectroscopy with multiple irradiation points for each multiplet using different expnos.
paropt	Parameter optimization au program.
parray	Parameter optimization au program using parameter arrays. Derived from 'paropt', but several parameters may now be changed per experiment. In addition, parameters are not changed via constant increments. Instead, the values are taken from an array.
pass2d	Perform a PASS experiment with 5 Pi-pulses and 16 increments (samples up to 16 spinning side bands).
pecosy	Program to pre-process P.E.COSY raw data before 2D-FT.
phtran	Transfer phase correction parameters PHC0 and PHC1 into acquisition parameters PH_ref and DE.
plintfac	Plot integrals with different scaling factors.
plot_sino	Plot spectrum, scaling depends on Signal/Noise.
plot_to_file	Creates a postscript file of the desired plot.
plotx	Plots individually scaled integral regions as separate objects.
poptau	Parameter optimization au program using parameter arrays. Derived from 'paropt' but several parameters can be optimized. The parameters are changed according to the parameter arrays. The AU program will be started from user interface 'popt' (parameter editor).
popthalt	Halt "popt" AU program.
proc_1H	Processes and plots 1D spectra. Does not perform baseline correction.
proc_1d	Processes and plots 1D spectra.

proc_1dapks	Processes and plot 1D spectra. Uses 'apks' for phase correction.
proc_1dconlf	Processes and plots 1D spectra. Plots an additional spectrum on the same plot if there are integrals in the lowfield range outside $\delta > 11$
proc_1dconlf_pr	Processes and plots 1D spectra. Plots an additional spectrum on the same plot if there are integrals in the lowfield range outside the plot limits.
proc_1dglp	Processing AU program with automatic data evaluation according to GLP standards. This AU program takes CY as an argument and then works like proc_1d.
proc_1dlf	Processes and plot 1D spectra. Plots an additional lowfield plot.
proc_1dpppti	Processes and plots 1D spectra. Creates a special peaklist file (frequency (Hz) and half width) and prints this on the plot.
proc_1dppti	Processes and plots 1D spectra. Creates a peak picking list and prints this on the plot.
proc_2d	Processing AU program for 2D spectra without plotting.
proc_2dhom	Processes and plots 2D homonuclear type spectra.
proc_2dhom_2pp	Processes and plots 2D homonuclear type spectra with two positive projections.
proc_2dinv	Processes and plots 2D inverse type spectra.
proc_2dinv_2p	Processes and plots 2D inverse type spectra. Plots two projections.
proc_2dphf2het	Determines phase correction in F2 for heteronuclear spectra.
proc_2dphf2hom	Determines phase correction in F2 for homonuclear spectra.
proc_2dpl	Processes and plots 2D type spectra.
proc_2dsym	Processes and symmetrizes 2D type spectra.

proc_2dt1	Automatic processing of one 2D T1/T2 experiment with subsequent T1/T2 calculation.
proc_cpd135	Processes and plots 13C CPD and DEPT135 spectra that were acquired with the AU program au_zg135.
proc_glp	Automatic GLP data evaluation.
proc_intrng	Processes and plots 1D spectra. Uses the predefined integral range file 'testrng' for integration.
proc_MAS	Processes and plot 1D MAS spectra.
proc_no	AU program which does no processing.
proc_noe	Processes and plots noediff spectra.
proc_t1	Semi-automatic processing of multiple 2D T1/T2 experiment with subsequent T1/T2 calculation.
proc_tecalib	Evaluation of previous temperature calibration experiments.
psys180f1t1	Processing AU program for the 180° pulse calibration tests.
psysamp1s39	Processing AU program for the amplitude stability tests - with shaped pulse - with 30 ° pulse - with 90° pulse - after gradient echo (5msec, 30 G/cm) - after gradient echo (5msec, 10 G/cm) - after gradient pulse (1msec, 10G/cm).
psysb1hom	Processing AU program for the B1 homogeneity test.
psysb2hom	Processing AU program for the B2 homogeneity test.
psyscancel	Processing AU program for the - phase cycling cancellation test - phase cycling cancellation test after gradient pulse.
psysdante1	Processing AU program for the dante type turn on test.
psysdecpro1	Processing AU program for the decoupler profile test.
psysexpro1	Processing AU program for the - excitation profile (16 usec gauss shape) test - excitation profile (6 msec gauss shape) test.
psysglitch	Processing AU program for the glitch test.

psysgreco1	Processing AU program for the gradient recovery test.
psysgrzpro	Processing AU program for the z-gradient profile.
psysmodl1	Processing AU program for the - modulator linearity test - shaped pulse modulator linearity test.
psysmulti1	Processing AU program for the amplitude linearity test (1dB power level steps).
psysphas1st	Processing AU program for the - phase stability test ("13° test") - shaped pulse phase stability test (16 usec gaussian shape, "13° test").
psysphasf1	Processing AU program for the - phase propagation test - phase shifting test.
psyspullin1	Processing AU program for the - amplitude linearity test - shaped pulse amplitude linearity test (pulse length *2, power level +6).
psysquadim	Processing AU program for the quad image suppression test.
psysrgtest	Processing AU program for the receiver gain test (analog and digital).
psyssoftp1	Processing AU program for the shaped pulse comparison (rectangular, gaussian, eburp1).

psystestab	Automatic processing of a 2D Temperature stability experiment, evaluation of temperature and statistic of temperature stability. Can be used to process data obtained with the AU program systestab.
psysturnon	Processing AU program for the turn on test.
pulse	Program to calculate attenuation value for given pulse length or nutation frequency, or vice versa.
pulsecal	Single scan pulse calibration via stroboscopic nutation experiment.
pulsecalib	Offset/Pulsecalibration H2O/D2O and Offset/Pulsecalibration 13C, 15N probehead.
qnpset	Define the QNP parameter according to the currently defined probehead.
quadplot	First plots a 2D overview spectrum and then the 4 quadrants of the 2D spectrum.
queue	Queue data acquisition.
queue_init	Initialise data acquisition with the AU program queue.
queuega	Queue data acquisition.
r23mplot	Read 2D slices from a 3D data set and plot them.
r23mult	Repeatedly reads slices from a 3D data set (3rrr) into successive experiment numbers.
rampXY	3D gradient shimming with the BSMS RCB board.
remproc	Automatic conversion and processing of data sets transferred via BRUKNET, LIGHTNET, NMR-LINK or TCP-LINK.
repeat	Repeat an acquisition with exactly the same parameters, pulse program and other lists.
rescale	Applies intensity scaling, for direct comparison of spectra acquire with different RG, NS and flip angle.
secplot	Generate a section plot. The overview spectrum is plotted together with a vertical expansion of a smaller part of the spectrum on top of it.
set2hdecgp	Setup AU program for standard 3D parameter sets.

setccnh3dgp	Setup AU program for standard 3D parameter sets.
setdiffparm	Extracts diffusion sequence parameters and stores parameters for "vargrad" simfit fitting (T1/T2) or DOSY processing.
seteditedgp	Setup AU program for standard 3D parameter sets.
sethccc3dgp	Setup AU program for standard 3D parameter sets.
setpar3dgp	Setup AU program for standard 3D parameter sets.
set_sreglist	Set SREGLST parameter from NUC1 and SOLVENT.
simplex	AU program for autoshimming. It is suitable for adjustment of strongly coupled shim groups which may be far from the optimum position.
simtoseq	Converts data which have been recorded in digital and qsim mode to data which appear to be acquired in qseq mode.
sinocal	Calculates the signal to noise ratio.
split	Separate data obtained with interleaved acquisition.
split2D	Splits a processed 2D file into single 1D spectra.
split3D	Separate interleaved 3D data.
splitcomb	Combine, shift and add 2D/3D data recorded with an interleaved single or double InPhase/AntiPhase or S3E scheme.
splitcrin	AU program to separate interleaved 3D data
splitcrinept	AU program to separate interleaved 3D data.
splitdqzq	Combine DQ/ZQ data.
splithb	Separate and combine H-bond experiments.
splithmsc	Separate the 1J and nJ component of the HMSC experiment.
splitinvnoe	Separate NOE and NONOE data obtained with a pulse program like invinoef3gpsi.
splitipap	Create separate InPhase and AntiPhase datasets.
splitipap2	Create separate InPhase and AntiPhase datasets for DSSE experiments.

splitpap3d	Separate Inphase/Antiphase data from 3D datasets.
splitpfids	Separate decoupled/non-decoupled data obtained with a pulse program like pfidsetgpsi.
splitser	Splits a ser file into single fids, starting with the expno which follows the ser file.
splitxf	Separate and combine double half filtered data.
ssd	Set the 2nd and 3rd dataset from the commandline.
stack1d	Generates a stacked plot of 1D spectra from increasing or decreasing EXPNOs or PROCNOs.
stack2d	Generate a 2D stack plot.
stackp1d	Generates a stacked plot of 3 to 12 1D spectra from increasing or decreasing EXPNOs or PROCNOs.
stdsplit	Splits STD pseudo 2D datasets.
suppcal	Calculates the width of the Water peak at 100% and 50% of the DSS signal height. The result is referred to as the 'water suppression test'.
sys180f1t1	Acquisition AU program for the 180° pulse calibration test with different phases.
sys180f1t2	Acquisition AU program for the 180° pulse calibration test with different flip angles.
sysamp1sp9	Acquisition AU program for the shaped pulse amplitude stability test.
sysamp1st	Acquisition AU program for the amplitude stability tests - with 30° pulse - with 90° pulse.
sysb1hom	Acquisition AU program for the B1 homogeneity test.
sysb2hom	Acquisition AU program for the B2 homogeneity test.
syscancel	Acquisition AU program for the phase cycling cancellation test.
sysdante1	Acquisition AU program for the dante type turn on test.
sysdecpro1	Acquisition AU program for the decoupler profile test.

sysexpro1	Acquisition AU program for the - excitation profile (16 usec gauss shape) test - excitation profile (6 msec gauss shape) test.
sysgenpar	Preparation AU program for all HWT test programs.
sysglitch	Acquisition AU program for the glitch test.
sysgrcan	Acquisition AU program for the phase cycling cancellation test after gradient pulse.
sysgrecho	Acquisition AU program for the amplitude stability test after gradient echo (5msec, 30 G/cm and 5msec, 10 G/cm).
sysgrreco1	Acquisition AU program for the gradient recovery test.
sysgrstab	Acquisition AU program for the amplitude stability test after gradient pulse (1msec, 10G/cm).
sysgrzpro	Acquisition AU program for the z-gradient profile.
sysmodl1	Acquisition AU program for the modulator linearity test
sysmodls1	Acquisition AU program for the shaped pulse modulator linearity test.
sysmultl1	Acquisition AU program for the amplitude linearity test (1dB power level steps).
sysphas1sp	Acquisition AU program for the shaped pulse phase stability test (16 usec gaussian shape, "13 degree test").
sysphas1st	Acquisition AU program for the phase stability test ("13 degree test").
sysphasf1	Acquisition AU program for the - phase propagation test - phase shifting test.
syspullin1	Acquisition AU program for the amplitude linearity test (pulse length *2, power level +6).
sysquadim	Acquisition AU program for the quad image suppression test.
sysrgtest	Acquisition AU program for the receiver gain test (analog and digital).

syssoftp1	Acquisition AU program for the shaped pulse comparison (rectangular, gaussian, eburp1).
sys spline1	Acquisition AU program for the shaped pulse amplitude linearity test (pulse length *2, power level +6).
systestab	AU program for a temperature stability experiment performed as pseudo 2D experiment including evaluation of temperature and statistics of temperature stability.
systemon	Acquisition AU program for the turn on test.
tecalib	AU program to determine the temperature calibration curve.
testsuite	Test the general functionality of a TOPSPIN release version. Basic functionality is given if this program is completed without error messages.
tmscal	Performs a peak picking around the TMS signal. If the two satellites from the ^{29}Si - ^1H coupling can be detected, the resolution is OK.
tune	Tune a probehead.
update_layout	Sets the parameter 'LAYOUT' in all parameter sets.
update_aunmp	Sets the parameter AUNMP in all parameter sets.
vtu_airflow	Switch the VTU AirFlow.
vtu_exam	Example AU program which shows how to use low level functions to read or write VTU (BVT1000/BDTC) parameters.
vtu_heater	Switch the VTU Heater.
writeshimZ	Reads the on-axis shim values and writes a pseudo shim file.
xfshear	Program for shearing of 2D MQMAS spectra of odd half integer quadrupolar nuclei. Data need to be acquired in States Mode
zeroim	Zeroe the imaginary data of a 1D or 2D data set.
zg_2Hoffon	General AU program for data acquisition. The lock is switched off before the acquisition is started.

zgchkte	Starts acquisition with zg and monitors the temperature. The experiment is halted if the current temperature differs too much from the target temperature.
zg_dfs	Calculates shape file for double frequency sweep and subsequent data-acquisition.
2df1shift	Shift a 2D spectrum along the F1 axis.
2dgetref	Gets parameters for a 2D spectrum from the 1D reference spectra : Nucleus, Frequencies, Spectral Width, and reference plot data set names. The F2 reference is taken from the second dataset. The F1 reference is taken from the third dataset.
2dshift	Shift 2D time domain data left or right over NSP points.
2nde	Set 2nd data set to new expno and 3rd data set equal to foreground data set.
2ndn	Set 2nd data set to new name and 3rd data set equal to foreground data set.

Chapter 14

TOPSPIN parameter types

This chapter contains a list of all TOPSPIN parameters grouped by their type. The type of a parameter can be integer, float, double or character-string. Several AU macros read TOPSPIN parameters into AU variables or store the value of AU variables into TOPSPIN parameters. In both cases it is important that the type of the AU variable is the same as the parameter type.

14.1 Integer parameters

The following TOPSPIN parameters are of the type integer:

ABSG	AQORDER	AQSEQ	AQ_mod
BC_mod	BYTORDA	BYTORDP	DATMOD
DIGMOD	DIGTYP	DS	EXPNO2
EXPNO3	FnMODE	FT_mod	HGAIN[4]
HOLDER	HPMOD[8]	HPPRGN	INTBC
L[32]	LOCSHFT	LPBIN	MC2
ME_mod	NBL	NC	NCOEF
NC_proc	NLEV	NS	NSP
NZP	OVERFLW	PARMODE	PH_mod
PKNL	POWMOD	PPARMOD	PRGAIN
PSCAL	PSIGN	PROCNO2	PROCNO3
QNP	REVERSE	RO	RSEL[10]
SI	STSI	STSR	SYMM
TD	TD0	TDeff	TDoff
TILT	WBST	WDW	XDIM
XGAIN[4]	YMAX_a	YMAX_p	YMIN_a
YMIN_p			

14.2 Float parameters

The following TOPSPIN parameters are of the type float:

ABSF1	ABSF2	ABSL	ALPHA
ASSFAC	ASSFACI	ASSFACX	ASSWID
AZFE	AZFW	BCFW	CNST[64]
DC	DE	D[64]	FCOR
FW	GAMMA	GB	GPX[32]
GPY[32]	GPZ[32]	INTSCL	ISEN
LB	LEV0	LOCPHAS	MASR
MAXI	MI	NOISF1	NOISF2
OFFSET	PC	PCPD[10]	PHC0
PHC1	PHCOR[32]	PH_ref	PL[64]
P[64]	RECPH	RG	SIGF1
SIGF2	SINO	SPOAL[32]	SPOFFS[32]
SP[32]	SSB	S_DEV	TE
TE2	TM1	TM2	TOPLEV
V9	VD		

14.3 Double parameters

The following TOPSPIN parameters are of the type double:

BF1	BF2	BF3	BF4
BF5	BF6	BF7	BF8
COROFFS	CY	F1P	F2P
INP[32]	IN[32]	LFILTER	LGAIN
LOCKPOW	LTIME	O1	O2
O3	O4	O5	O6
O7	O8	SF	SFO1
SFO2	SFO3	SFO4	SFO5
SFO6	SFO7	SFO8	SW
WBSW			

14.4 Character-string parameters

The following TOPSPIN parameters are of the type character-string:

AUNM[16]	AUNMP[16]	CPDPRG[16]	DFILT[16]
DSLST[16]	DU[256]	DU2[256]	DU3[256]
EXP[32]	FQ1LIST[16]	GPNAM0[64]	GRD- PROG[16]
INSTRUM[64]	LAYOUT[256]	LOCNUC[8]	MASRLST[16]
NAME[64]	NUC1[8]	PROBHD[64]	PULPROG[32]
SOLVENT[32]	SPNAM0[64]	SREGLST[40]	TI[72]
TYPE[16]	USER[64]	USERA1[80]	USERP1[80]
VCLIST[16]	VDLIST[16]	VPLIST[16]	VTLIST[16]

Index

A

ABORT 46
ABS 33
ABS1 37
ABS2 37
ABSD 33
ABSD1 37
ABSD2 37
ABSF 33
ABSF1 33
ABSF2 33
ABSOT1 37
ABSOT2 37
ABST1 37
ABST2 37
ADD 35, 63
ADD_CURDAT_TO_PORTFOLIO 43, 114
ADD_TO_PORTFOLIO 43, 115
ADD2D 37
ADDC 35
addfid command 104
AND 35
APK 33, 70
APK0 33
APK1 33
APKF 33
APKS 33
AU macro 6
aucmd.h 16, 19
AuditAppend 130, 131
AUDITCOMMENTA 26
AUTOGAIN 29
Automatic baseline correction 41
AUTOPHASE 29
AUTO PLOT 42, 110
AUTO PLOT_TO_FILE 42, 112
AUTO PLOT_WITH_PORTFOLIO 43, 117
AUTO PLOT_WITH_PORTFOLIO_TO_FILE 43,

118

AUTOSHIM_OFF 29
AUTOSHIM_ON 29
autoshimming 30

B

base_info file 35
baseline correction
 ID user defined 33
baslpnts file 35
BC 33
BCM1 37
BCM2 37
brukdef.h 19
Bruker library functions 9

C

CalcExpTime 126
CalcExpTime function 60
cc compiler 16
C-code 16
character string parameters 178
ChecksumFile 127
CLOSE_PORTFOLIO 43, 116
column of a 2D spectrum 79, 82
compileall command 7, 8
compiling AU programs 7
constants 16
control statements 15
CONVDTA 33
cplbruk command 8
CPR_exec 15, 44, 50, 51, 60
CREATE_PORTFOLIO 43, 113

D

DATASET 25, 60, 62, 71
DATASET2 25, 63
DATASET3 25, 63
DDATASETLIST 25

DECLARE_PORTFOLIO 43
define statements 13, 16
DEG90 29
DELPAR 28
DEXPNO 25, 61, 66
dircp 137
dircp_err function 137
DIV 35
dpa command 98
DPROCNO 25, 69
DPULPROGLIST 29
DT 35
DU 26
DVTLIST 31

E

eda command 44, 96, 98, 106, 107
edau command 7, 8, 10
edc2 command 25, 115
eddosy command 28
edit mode 8
edlock command 30
edmac command 9
edmisc command 37
edo command 106, 107
edp command 44, 96, 106
EF 33
EFP 33, 64
EJ 29
EM 33
enhanced metafile 112, 118
erropt.h 19
ERRORABORT 46
Executing AU programs 8
expinstall command 7

F

F1DISCO 40
F1PROJN 40
F1PROJP 40
F1SUM 40
F2DISCO 40
F2PROJN 40
F2PROJP 40
F2SUM 40

fcntl.h 19
FETCHDOSYPAR 28
FETCHPAR 27, 96
FETCHPAR1 27, 96
FETCHPAR1S 27
FETCHPAR3 27, 96
FETCHPAR3S 27
FETCHPARM 28, 96
FETCHPARN 27
FETCHPARNS 27
FETCHPARS 27, 98
FETCHPARS1 98
FETCHPARS3 98
FETCHT1PAR 28, 96
fidtoser AU program 86
FileSelect 132
FILT 35
first order phase correction 33
float parameters 177
FMC 33
FP 33
freedir 136
FROMJDX 44, 92
FT 33, 70
ft command 9

G

Gaussian deconvolution 37
Gaussian window multiplication 34
gcc compiler 16
GDATASETLIST 25
GDCHECK 26
GDCHECKRAW 26
GDCON 37
GENFID 34
GENSER 39
GETCURDATA 20, 25
GETCURDATA2 25
GETCURDATA3 25
GETDATASET 25, 64
getdir 134
GETDOUBLE 27, 121
GETFLOAT 27, 99, 121
gethighest 138
GETINT 27, 120

GETLCOSY 42
GETLIM 42
GETLINV 42
GETLJRES 42
GETLXHCO 42
getParfileDirForRead 139
getParfileDirForWrite 142
GETPROSOL 27
getstan 144
GETSTRING 27, 123
getxwinvers 146
GF 34
GFP 34
GLIST 25
GM 34
GO 29
GPULPROGLIST 29
GVTLIST 31

H

header files 19
Hilbert Transform 34, 39
HT 34

I

IDATASETLIST 21, 25
IEXPNO 25, 51, 60, 65, 71
IFEODATASETLIST 25
IFEOPULPROGLIST 29
IFT 34
II 29
IJ 29
ILOOPCOUNTLIST 21
include statements 13, 14, 16
integer parameters 176
intrng file 35, 123
Inverse Fourier Transform 34
 2D 39
INVSF 38
IPARSETLIST 21
IPROCNO 25, 68
IPULPROGLIST 21, 29
IVTLIST 21, 31

J

JCAMP-DX file 44, 92
JCAMP-DX format 35, 90
JCONV 44, 94
Jeol dataset 44, 94

L

lastparflag variable 11
LDCON 37
LEVCALC 38
LFILTER 30
LG 30
LGAIN 30
LI 35, 123
libcb.h 19
limits.h 19
LIPP 35
LIPPF 35
listall_au AU program 9
LO 30
LOCK 30
lock power 30
LOCK_OFF 30
LOCK_ON 30
LOCNUC 30
loop gain 30
loop statements 15
loop structures 13
loop time 30
loopcount1 variable 11
loopcount2 variable 11
LOPO 30
LS 35
LTIME 30

M

Magnitude calculation 34
MAKE_ZERO_FID 29
makeau file 16
MAS unit 31
MASE 31
MASG 32
MASH 32
MASI 31
MASR 31

MASRGET 32
math.h 19
MC 34
MDCON 37
mkudir 147
MUL 36, 63
MULC 36
multiexpt AU program 126
multizg AU program 60

N

NM 36
NZP 36

P

parameter type 175
PathXWinNMR function 148
peaklist file 35
phase correction first order 33
phase correction zero order 33
PHC0 33
PHC1 33
PK 34
plane from 3D raw data 87
plot_to_file AU program 110
portfolio of XWIN-PLOT 113, 114, 115, 116, 117,
118
postscript file 110, 112, 118
pow_next 150
Power spectrum 38
PP 35
PPH 35
PPP 35
ppp command 37
predefined dedicated variables 10
predefined general variables 10
PrintExpTime 126
PrintExpTime function 60
Proc_err 151
Proc_err function 14
processed data 74, 75, 76
PS 34
PTILT 38
PTILT1 38

Q

QSIN 34
quick reference 7
QUIT 46
QUITMSG 46

R

R12 41
R13 41
R23 41
raw data 74, 76
RDATASETLIST 25
reg file 35
REV1 38
REV2 38
REXPNO 25, 67
RGA 29
RHNP 40
RHPP 40
RMISC 35, 123
ROT 30
rotation 30
ROTOFF 30
row of 2D raw data 83, 85
row of a 2D spectrum 78, 81
RPAR 28, 106
RPROCNO 25, 70
RPULPROGLIST 20, 29
RS 36
RSC 40, 79
RSER 41, 83
RSER2D 87
RSH 30
RSR 40, 78
RV 36
RVNP 40
RVPP 40
RVTLIST 20, 31

S

SAB 34
sample.h 19
second AU dataset 63
SETCURDATA 25, 60, 61
SETDATASET 26

SETPULPROG 29
SETSH 30
SETUSER 26
Show_status 153
showfile 154
Sine window multiplication 34
SINM 34
SINO 34
SOLVENT 30
Spline baseline correction 34
splitser AU program 84
SREF 34
ssleep 155
stdio.h 17, 19
stdlib.h 17, 19
STOP 46, 47
STOPMSG 46, 47
STOREDOSYPAR 28
STOREPAR 27, 100
STOREPAR1 27, 100
STOREPAR1S 27, 104
STOREPAR3 27, 100
STOREPAR3S 27, 104
STOREPARM 28
STOREPARN 27, 102
STOREPARNS 28
STOREPARS 27, 104
STORET1PAR 28, 100
strcpy C-function 62
string.h 19
SUB1 38
SUB1D1 38
SUB1D2 38
SUB2 38
subroutines 10, 11
SWEEP_OFF 30
SWEEP_ON 30
SYM 38
SYMA 38
SYMJ 38

T

TABS1 41
TABS2 41
TABS3 41

Tcl/Tk scripts 6
TE2GET 31
TE2READY 31
TE2SET 31
TEGET 31
temperature unit 31
TEPAR 31
TEREADY 31
TESET 31
TF1 41
TF1P 41
TF2 41
TF2P 41
TF3 41
TF3P 41
third AU dataset 63
TILT 38
TIMES2 19
TIMES3 19
TIMESLIST 21
TM 34
TOJDX 44, 90
TOJDX5 90
Trapezoidal baseline correction 37
Trapezoidal window multiplication 34
TRF 34
TUNE 30
TUNESX 30

U

uni.h 19
unistd.h 19
unlinkpr 156
USECURPARS 19
USELASTPARS 19
user defined
 baseline correction 33
user defined variables 10, 11
util.h 19
UWM 34

V

variable assignments 13
variable declarations 13
Varian dataset 44, 93

VCONV 44, 93
view mode 8
VIEWDATA 26, 66, 69, 71, 72
viewing AU programs 8
vorspann file 16
VT 31
VTLIST 31

W

WAIT_UNTIL 44, 55
WMISC 35
WPAR 28, 107
WRA 26, 74
WRP 26, 75
WRPA 26, 76
WSC 41, 82
WSER 41, 85
WSERP 41
WSH 30
WSR 81

X

XAU 44
XAUA 44
XAUP 44
XAUPW 44
XCMD 44, 54
XF1 38
XF1M 38
XF1P 38
XF1PS 38
XF2 38
XF2M 38
XF2P 38
XF2PS 38
XFB 38
XFBM 38
XFBP 38
XFBPS 38
XHT1 39
XHT2 39
XIF1 39
XIF2 39
XMAC 44
xmac command 9

XTRF 39
XTRF2 39
XTRFP 39
XTRFP1 39
XTRFP2 39
XWP_LP 42
XWP_PP 42

Z

zero order phase correction 33
ZERT1 39
ZERT2 39
ZF 36
ZG 29
zg command 9
ZP 36